SURGICEL™

ENDOSCOPIC APPLICATOR

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






SURGICEL™ ENDOSCOPIC APPLICATOR

en

 Do not use device intravascularly because as with any topical haemostatic agents introduced into the circulatory system, life-threatening embolic events may occur.

 To prevent clogging, do not touch the tip to wet surface. Be careful to avoid damaging tissue with the rigid tip.

READ INSTRUCTIONS FOR SURGICEL™ POWDER PRIOR TO USE

DESCRIPTION

The SURGICEL™ Endoscopic Applicator is supplied with a flexible inner tip and a rigid cannula. The SURGICEL™ Endoscopic Applicator is disposable and must be discarded after use.

INDICATIONS

The SURGICEL™ Endoscopic Applicator is intended for use in delivering powdered haemostatic agent to bleeding surgical sites through a 5 mm or larger trocar.

WARNINGS

- To prevent clogging, do not touch the tip to wet surface. Be careful to avoid damaging tissue with the rigid tip.
- Do not attempt to trim the applicator tip. Replace the tip if it becomes clogged.
- Do not use the SURGICEL™ Endoscopic Applicator if package is opened or damaged.
- Apply SURGICEL™ Powder according to the product's labeling.
- Do not resterilize/reuse. Reuse of this device (or portions of this device) may create a risk of product degradation, which may result in device failure and/or cross-contamination, which may lead to infection or transmission of blood-borne pathogens to patients and users.

PRECAUTIONS

- The SURGICEL™ Endoscopic Applicator is supplied with a flexible inner tip inside a rigid cannula. The rigid cannula cannot be used independently.
- The SURGICEL™ Endoscopic Applicator should only be used by persons having adequate training and familiarity with endoscopic techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any endoscopic procedure.
- After use, dispose of the SURGICEL™ Endoscopic Applicator in accordance with biohazardous material protocol.

DIRECTIONS FOR USE

1. Open package and transfer SURGICEL™ Endoscopic Applicator to sterile field. The SURGICEL™ Endoscopic Applicator is composed of two parts, a flexible inner tip and a rigid cannula.
2. Remove the tip from the SURGICEL™ Powder Device (Fig. 1).
3. Twist to open (Figs. 2, 3).
4. Attach the SURGICEL™ Endoscopic Applicator to the SURGICEL™ Powder Device without pumping the bellows and ensure the connection is secure (Fig. 4).
5. The device with the rigid tip is now ready for use (Fig. 5).
6. If a completely flexible tip is desired, remove the rigid cannula (Fig. 6 - optional).
7. Introduce the SURGICEL™ Endoscopic Applicator into the trocar port.

CAUTION: Insert the SURGICEL™ Endoscopic Applicator using caution not to inadvertently express powder.

THIS PAGE IS LEFT BLANK INTENTIONALLY.

8. Carefully position the distal end of the SURGICEL™ Endoscopic Applicator at the site where the powder is to be delivered. To prevent clogging, do not touch the tip to wet surface. Be careful to avoid damaging tissue with the rigid tip.
9. While holding the SURGICEL™ Endoscopic Applicator in place, pump to apply powder to treatment site.
10. Carefully remove the SURGICEL™ Endoscopic Applicator from the trocar port when sufficient powder has been delivered to the treatment site. Apply direct pressure over the treated site as needed.

STORAGE

No special storage conditions required. Do not use after expiration date.

HOW SUPPLIED

The SURGICEL™ Endoscopic Applicator is supplied sterile in a sealed package. Sterility of contents guaranteed unless package is opened or damaged.

SYMBOLS USED ON LABELING

| | |
|--|---|
|  Do not reuse |  Caution |
|  Do not use if package is damaged |  Catalogue number |
|  Do not resterilize |  Sterilized using Irradiation |
|  Use-by date |  Authorized Representative in the European Community |
|  Manufacturer |  CE Mark and Identification Number of Notified Body. The product meets the essential requirements of Medical Device Directive 93/42/EEC. |
|  Batch code |  Quantity |

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



IFU PRINTING SPECIFICATION SHEET

| | | | | | | | | | | |
|---|--|---|---|---|--|--|---|--|------------------------------|---|
| TITLE SURGICEL™ Endoscopic Applicator CE | | DESCRIPTION Booklet IFU with flap | | LAB NUMBER LAB100481641v3 | SPECIAL INSTRUCTIONS/COMMENTS 1 panel fold in from front cover, Staple text pages into cover as shown | BINDING Saddle Stitch | COLORS Black | | | |
| FLAT SIZE (COVER / PAGES) 12.5" x 9.625" (317.5 mm x 244.475 mm) 10" x 9.625" (254 mm x 244.475mm) | FOLDED SIZE (COVER / PAGES) 5" x 9.625" (127 mm x 244.475 mm) 5" x 9.625" (127 mm x 244.475 mm) | RMC NUMBER 389960R03 | PAGE COUNT 60 + 6pg cover | LANGUAGES ar, cs, da, de, el, en, es, et, fi, fr, hr, hu, it, ko, lt, lv, nl, no, pl, pt, ro, ru, sk, sl, sr, sv, tr, zh-cn, zh-tw | | | SELF COVER <input type="checkbox"/> | PLUS COVER <input checked="" type="checkbox"/> | SEALING METHOD n/a | WAFER SEAL <input type="checkbox"/> |
| BLEED SIZE .5" (12.7 mm) <input type="checkbox"/> .125" (3.175 mm) <input type="checkbox"/> | NONE <input checked="" type="checkbox"/> BLEED ALL SIDES <input type="checkbox"/> | BLEED TOP <input type="checkbox"/> | BLEED RIGHT <input type="checkbox"/> | BLEED LEFT <input type="checkbox"/> | BLEED BOTTOM <input type="checkbox"/> | DRAWING IS NOT TO SCALE: DRAWINGS REFLECT INFORMATION FOR PRODUCTION OF PRINTED PIECES AND DO NOT CONTAIN ACTUAL ARTWORK. This document or data herein or herewith is not to be reproduced, used or disclosed in whole or part without the permission of Ethicon, Inc. | | | | |
| STOCK 50# cover finch opaque, 30# text | | | | | ETHICON | | | | | |

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© Ethicon, Inc. 2017