# MR370

MR700  
MR730  
MR810  
MR850

MR410

English

(en)

## Introduction

The MR370 is a reusable humidification chamber for use with Fisher & Paykel heater bases and is recommended for adult applications. The chamber is used to provide humidification for intubated or non-intubated patients.

## Product Specifications

Gas Ports:	Inlet	22 mm Female
	Outlet	22 mm Male
Maximum Operating Pressure:	20 kPa	
Compressible Volume: Full	With scroll	430 mL
	Without scroll	430 mL
Empty	With scroll	700 mL
	Without scroll	715 mL
Maximum Water Capacity:	280 mL	
Compliance:	0.6-0.8 mL/cm H <sub>2</sub> O	
Maximum Peak Flow:	180 L/min	
Materials Used In Construction:	Base	Aluminium
	Dome	Polysulphone
	Gasket	Synthetic Rubber

## Water level symbol

- x denotes water level too high.
- ✓ denotes maximum water level.

## Set Up and Operation

- Ensure that the chamber and circuit are sterilised before use.
- If an MR700 series, MR730 or MR850 humidifier is used, the scroll is not required. If an MR410 or older series humidifier is used, load the paper into the scroll of the chamber following the instructions included with the paper.
- Refer to the humidifier operating manual for recommended set-up on the humidifier.
- Replace and sterilise the chamber in accordance with hospital infection control procedures.

## Warning

- DO NOT fill the chamber above the maximum fill level line. Liquid could enter the breathing circuit if the chamber is overfilled.
- DO NOT fill the chamber with water in excess of 37 °C.
- If waterfed set is used, ensure the water bag is placed below the level of the humidifier after the chamber has been filled.
- Ensure that the water level in the chamber is periodically monitored. Refill when necessary.
- Ensure the humidifier is not tilted. Tilting the humidifier may result in water entering the breathing circuit.
- DO NOT touch the heater plate or chamber base. Surfaces may exceed 85°C.
- Do not exceed a maximum peak flow rate of 180 L/min.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: [www.fphcare.com/prop65](http://www.fphcare.com/prop65)

## Cleaning

Chambers may be cleaned by any of the following methods:

Autoclave:	132°C (270°F) @ 220kPa (32psi) for 4 minutes
	or 121°C (250°F) @ 96kPa (14.1psi) for 30 minutes
Ethylene Oxide:	55°C (131°F)

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Chemical Immersion: Solutions such as Medizyme, Pyroneg, Control 3, Solution 2 and Cidex have been tested according to the manufacturers directions. Other brand names with similar active ingredients may also be suitable.

## Warning

The following solutions should be avoided as they may cause the chamber to crack:

- Cetones, Formaldehyde, Chlorinated Hydrocarbons, Hypochlorite, Inorganic Acids, Aromatic Hydrocarbons, Phenol (>5%)

To prevent cracking, ensure that the water inlet port plug, and any other reusable components, are disconnected from the chamber before cleaning.

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REF 185041322 REV L 2019-10 © 2019 Fisher & Paykel Healthcare Limited