of the quality control material *EasyCon_{norm}*

Intended Use:

EasyCon is a liquid quality control material intended for use in general chemistry to monitor the precision and accuracy of glucose- and lactate measuring.

Handling Instructions:

EasyCon is recommended to be used on the following devices: BIOSEN C_line and BIOSEN S_line.

The quality control material **EasyCon** must be prepared and evaluated in the same way as a patient sample. Which means that a 20µL end-to-end capillary must be filled with the material (the capillary can be filled directly out of the storage container) and given after that into a pre-filled reaction cup, shook and placed in the device to be evaluated. After use of the storage container it should be closed properly and placed back in the fridge.


Warning:

Before applying, please refer to devices manual chapter „Preparing samples„. The quality control material must be diluted before evaluating in BIOSEN devices (see Handling).

Treat the material as **potentially infectious**.

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Linearity test kit, glucose / lactate, 2 mmol/L, 7 mmol/L, 18 mmol/L
Diluted 51-fold, ready for use

Intended Use:

For in vitro diagnostic use with BIOSEN C_line and BIOSEN S_line glucose and lactate analyzers. Also for non-diagnostic use. Not for food or drug use.

Handling Instructions:

Set-Lin-Test Glucose/Lactate can be used for linearity or calibration verification studies. Place the micro test tube in a sample position of the sample tray and start measurement as usual.

Chemical Composition:

0,039 mmol/L Lactate (2 mmol/L diluted 51-fold),
 0,137 mmol/L Lactate (7 mmol/L diluted 51-fold) and
 0,353 mmol/L Lactate (18 mmol/L diluted 51-fold)
 in a phosphate buffered solution (pH 7.2) containing further ingredients.

Warning:

Avoid contact with eyes and skin. Do not swallow.
 Keep away from foodstuff, beverages and food.

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Values:

REF	Item-No.:	0209-0102-391	3x 1 mL 2 mmol/L in micro test tubes blue 3x 1 mL 7 mmol/L in micro test tubes green 3x 1 mL 18 mmol/L in micro test tubes red
LOT	Charge-No.:		see label
	Expiry date unopened		see label
			Discard if expiration date is past.
	Lagertemperatur:	2 – 8 °C	
	Stability after opening:	8 hours (15 – 25 °C)	

Mean Glucose	Range
2.00 mmol/L 36.0 mg/dL	1.80 – 2.20 mmol/L 32.4 – 39.6 mg/dL
7.00 mmol/L 126.1 mg/dL	6.75 – 7.25 mmol/L 121.6 – 130.6 mg/dL
18.00 mmol/L 324.1 mg/dL	17.50 – 18.5 mmol/L 315.1 – 333.1 mg/dL

Mean Lactate	Range
2.00 mmol/L 18.0 mg/dL	1.80 – 2.20 mmol/L 16.2 – 19.8 mg/dL
7.00 mmol/L 63.1 mg/dL	6.75 – 7.25 mmol/L 60.8 – 65.3 mg/dL
18.00 mmol/L 162.1 mg/dL	17.50 – 18.5 mmol/L 157.6 – 166.7 mg/dL

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Intended Use:

In-vitro-diagnostics

Sensor test solution, glucose and lactate, ready for use, is used to check the selectivity of chip sensor glucose type II and chip sensor lactate type II.

Handling Instructions:

Sensor test solution is recommended to be used on the following devices: BIOSEN C_line and BIOSEN S_line.

You should measure **sensor test solution** daily as a control sample on the first 3 days of sensor use and after that at least twice a week. The results obtained with that solution must be within 4.8 mmol/L (86 mg/dL glucose, 43 mg/dL lactate) and 6.4 mmol/L (115 mg/dL glucose, 57 mg/dL lactate). Otherwise the sensor must be replaced. That solution detects the mechanical integrity of the sensor.

Chemical Composition:

5.0 mmol/L glucose, 5.0 mmol/L lactate, 1 mmol/L uric acid and 1 mmol/L acetaminophen diluted 51-fold a phosphate buffered solution (pH 7.2) containing further ingredients.

Warning:

Avoid contact with eyes and skin. Do not swallow.

Keep away from foodstuff, beverages and food.

After skin contact rinse with water. After eye contact rinse opened eye under running water. If symptoms persist, consult doctor.

Details:

	Item No.:	5130-6055	Sensor test solution, glucose and lactate, ready for use, set of 20 pcs.
			see label
	Expiry date unopened tubes:		see label
	Storage temperature:	2 °C - 8 °C	

Instructions for use

Multi Standard Solution 12 mmol/L (216.2 mg/dL) glucose and
12 mmol/L (108.7 mg/dL) lactate
diluted 51-fold , ready for use

Intended Use:

For in vitro diagnostic use with BIOSEN, EBIO, ESAT and ECONOMIC enzymatic instruments.

Also for non-diagnostic use.

Not for food or drug use.

Handling Instructions:

See instruction booklet of the instrument.

Chemical Composition:

0.235 mmol/L glucose (12 mmol/L diluted 51-fold) and 0.235 mmol/L lactate (12 mmol/L diluted 51-fold) in a phosphate buffered solution (pH 7.2) containing further ingredients.

Traceability:

The value of the calibrator for glucose is traceable to NIST (National Institute of Standards & Technology, Gaithersburg USA) Standard Reference Material 917b (D-Glucose).

Concentration: 12 mmol/L \pm 0.25 mmol/L, diluted 51-fold, ready for use.

Warning:

Avoid contact with eyes and skin. Do not swallow.

Keep away from foodstuff, beverages and food.

After skin contact rinse with water. After eye contact rinse opened eye under running water. If symptoms persist, consult doctor.

Values:

to REF (Item No.):	5211-3011	10x 2mL in micro test tube red
to REF (Item No.):	5211-3015	50x 2mL in micro test tube red
to REF (Item No.):	5211-3017	100x 2mL in micro test tube red
to REF (Item No.):	5211-3028	50mL in bottle
LOT:	see label	
⌚ (Expiry date unopened):	see label	
Storage:	5 – 8 °C Discard if expiration date is past.	
Stability after opening:	8 hours (15 –25°C)	

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Test serum solution ReadyCon_{Pat} for glucose and lactate

Intended Use:

Liquid, pre-diluted 1:51-fold, ready for use test serum with glucose and lactate concentrations in the pathological range.

For in vitro diagnostic use with the following glucose and / or lactate analyzers: BIOSEN C_line and BIOSEN S_line. Use the product only in combination with reagents from EKF-diagnostic GmbH.

Attention: simultaneous use of reagent from different suppliers on one instrument can lead to erroneous results.

Also for non-diagnostic use. Not for food or drug use

Handling Instructions:

The ready for use test solution must be placed in the sample or control solution positions and must be measured separate or together with the patient samples.

Warning:

Please note that at use of a ready for use pre-diluted test solution the filling of the capillaries and the quality and volume of the hemolyzing solution is not controlled.

Before applying, please refer to operator's manual of the device chapter „Measurement operation“. Treat the material as **potentially infectious**.

Values:

REF	Item No.:	5130-6162	25x 1mL in 2.0 mL brown safe lock micro test tubes
LOT		XXX	
	Expiry date unopened:	JJJJ-MM	 8°C
	Stability after opening:	8 hours (15 – 25°C)	2°C 
	Target value glucose:	16.38 mmol/L	295.2 mg/dL
	Confidence limits:	14.58-18.18 mmol/L	262.7-327.7 mg/dL
	Target value lactate:	15.29 mmol/L	137.6 mg/dL
	Confidence limits:	13.61-16.97 mmol/L	122.5-152.7 mg/dL

8001-0121-2300_gaEY__Testserum-ReadyCon-Pat_Gebrauchsanweisung_en+ru_A4-2xA5_20JJ-MM-LOT__VORLAGE