# ABLUMINUS DES+

## SIROLIMUS ELUTING CORONARY STENT SYSTEM

EN	Sirolimus eluting coronary stent system Instructions for use	Page 04
FR	Système de stent coronaire à élution de sirolimus Instructions d'utilisation	Page 09
DE	Sirolimus-freisetzendes Koronares stent system Gebrauchsanweisung	Page 14
IT	Sistema di stent coronarico a eluizione di Sirolimus Istruzioni per l'uso	Page 20
ES	Sistema de Stent CoronarioLiberador de Sirolimus Instruccionesdeuso	Page 26
PT	Sistema de Stent Coronário com Eluição de Sirolimus Instruções de uso	Page 32
RU	Сиролимус-элютирующая коронарная стентная система Инструкции по использованию	Page 37
TR	Sirolimus Salınlımlı Koroner Stent Sistemi Kullanım talimatları	Page 43
SR	Sirolimus Eluting Coronari Stent Sistem Uputstva za upotrebu	Page 49

## Symbol used in labeling

Length  
Longueur  
Länge  
Lunghezza  
Longitud  
Comprimento  
длина  
Uzunluk  
Dužina

Do not re-use  
Ne pas réutiliser  
Nicht wiederverwenden  
Non riutilizzare  
Não reutilizar  
Tekrar kullanmayın  
No reutilizar  
Не использовать повторно  
Nemojte ponovo koristiti

Diameter  
Diamètre  
Durchmesser  
Diametro  
Diámetro  
Diámetro  
диаметр  
Çap  
Prečnik

Use by Date  
Date limite d'utilisation  
Herstellungsland  
Data di scadenza  
fecha de caducidad  
Data de validade  
срок годности  
Son Kullanma Tarihi  
Koristiti po datumu

### REF

Catalogue Number  
Numéro de catalogue  
Katalognummer  
Numero del catalogo  
Número de catálogo  
Número do Catálogo  
Каталожный номер  
Katalog Numarası  
Kataloški broj

### STERILE EO

Sterilized using Ethylene Oxide  
Stérilisé à l'oxyde d'éthylène  
Sterilisiert mit Ethylenoxid  
Sterilizzato con ossido di etilene  
Esterilizado con óxido de etileno (EtO)  
Esterilizado com óxido de etileno  
Стерилизовано с помощью этиленоксида  
Etilen Oksit ile sterilize edilir  
Sterilizovano korišćenjem etilen oksida

### LOT

Batch Code  
Code du lot  
Code der Charge  
Codice del lotto  
Codigo de lote  
Código de lote  
Код серии  
Seri kodu  
Serijski kod

Caution,  
Attention  
Vorsicht,  
Attenzione,  
¡Precaución!  
Precaução  
Внимание,  
Dikkat,  
Oppez

### SN

Serial Number  
Numéro de série  
Seriennummer  
Numero di serie  
Número de serie  
Número de Série  
серийный номер  
Seri Numarası  
Serijski broj

### MD

Medical device  
Dispositif medical  
Medizinisches Gerät  
Dispositivo medico  
Dispositivo médico  
Dispositivo Médico  
Медицинское устройство  
Tıbbi cihaz  
Medicinski aparat

## Symbol used in labeling

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**1.0 PRODUCT DESCRIPTION**

The Abluminus DES+ Sirolimus Eluting Coronary Stent System is a device/drug combination product consisting of two major components: a) a bare metal stent as an implant (L-605 Cobalt Chromium Alloy coronary stent) crimped on a standard rapid exchange delivery system, and b) a unique coating that is applied in two layers on the crimped stent system (stent/balloon), with the first active layer consisting of a formulation containing Sirolimus and a second layer consisting of Polyvinyl Pyrrolidone as a protective top coat without Sirolimus. Generic name of Abluminus DES+ is "Sirolimus Eluting Coronary Stent System" with GMDN code 58771.

**1.1 DEVICE COMPONENT DESCRIPTION**


Abluminus DES+ is a balloon expandable intracoronary stent pre-mounted onto a rapid exchange delivery catheter. Three main components are (a) bare metal stent, (b) rapid exchange delivery system and (c) drug and polymer matrix coating. Each component is described below.

**(a) Bare metal Stent**


The bare metal stent of Abluminus DES+ is L-605 Cobalt Chromium Alloy. Stent design is sinusoidal two connector hybrid stent design with 6 cells, 8 cells and 10 cells design for in 2.25 mm - 2.50 mm, 2.75 - 3.50 mm and 4.00 - 5.00 mm respectively. The stent hybrid stent design provides excellent radial strength, flexibility and low crossing profile with optimal deployment characteristics for homogeneous expansion.

**(b) Rapid exchange delivery system**

The bare metal stent is crimped onto a rapid exchange semi-compliant balloon catheter located near distal tip. The distal shaft of the catheter has two lumens, one for inflation of the balloon and another for a guide wire advancement through the stenosis to be stented. The balloon part of catheter is made of polyamide and provides an expandable segment of a known diameter at a nominal pressure. The proximal shaft is made from a stainless steel hypotube. Proximal visual markers located 90 cm and 100 cm from the distal tip to aid catheter positioning using either the brachial or femoral approach. The catheter working length is 140 cm and is available in diameters between diameter of 2.25 mm and 5.00 mm, and stent mounted on balloon lengths ranging from 8 mm to 52 mm, as shown in Table 1. above. A schematic drawing of the Abluminus DES+ is shown in Figure 1.

**Figure 1: Schematic of Abluminus DES+ Sirolimus Eluting Coronary Stent System**

**Figure 2: Cell design of Schematic of Abluminus DES+**

**Table 1: Device Component Description**


Stent Design	6 cells	8 cells	10 cells
Stent Diameter (mm)	2.25-2.50	2.75-3.50	4.00-5.00
Stent length (mm)	8, 12, 16, 20, 24, 28, 32, 36, 40, 44, 48, 52		
Stent Material	L-605 alloy (Cobalt Chromium)		
Nominal Pressure	8 Bar		
Rated Burst Pressure*	14 Bar		
Delivery System Workable Length	140cm (1400mm)		
Catheter Configuration	Rapid Exchange		
Guidewire compatibility (max.)	0.014"		
Guiding Catheter Inner Diameter	5F (1.65mm) (inner lumen ≥ 0.058")		
Radiopaque Marker Bands	90% Pt / 10% Ir alloy (2)		

Formulation	Sirolimus drug, an active pharmaceutical ingredient in a bio-degradable polymer matrix
Sterilization	Ethylene Oxide Gas

\*Do not exceed Rated Burst Pressure

(c) Drug and Polymer matrix coating C1) DRUG: The active pharmaceutical ingredient on the Abluminus DES+ is Sirolimus. Sirolimus is a white powder, produced by a Streptomyces Hygroscopicus, with the empirical formula  $C_{51}H_{85}NO_{13}$  and a molecular weight of 914.2g/mole. It is insoluble in water and melts at approximately 183°C -185°C. The chemical name for Sirolimus is: (3S, 6R, 7E, 9R, 10R, 12R, 14S, 15E, 17E, 19E, 21S, 23S, 26R, 27R, 34aS)-9, 10, 12, 13, 14, 21, 22, 23, 24, 25, 26, 27, 32, 33, 34, 34a - Hexadecahydro - 9, 27 - dihydroxy - 3 - [(1R) - 2 - [(1S, 3R, 4R) - 4 - hydroxy - 3 - methoxycyclohexyl] - 1 - methylethyl] - 10, 21 - dimethoxy - 6, 8, 12, 14, 20, 26 - hexamethyl - 23, 27 - epoxy - 3H - pyrido [2, 1 - c] [1, 4] oxazacycloheptanone - 1, 5, 11, 28, 29 (4H, 6H, 31H)-pentone; Sirolimus CAS Registry number is 53123-88-9.

Sirolimus has the following chemical structure:

**C2: Biodegradable Polymers:**

Lactide and Glycolide family of bio-degradable polymers are used as drug carrier to control drug release from the stent.

**Formulation:**

A combination of bio-degradable polymers mixed with sirolimus makes up the coating formulation which is applied on the abluminus surface of Stent and exposed parts of balloon in pre-crimped configuration.

Nominal dosages of Sirolimus on Abluminus DES+ Stent with stent sizes is mentioned in table 2.

**Table 2: Product Matrix and Sirolimus Content**

Product Reference	Nominal Expanded Stent Diameter (mm)	Nominal Expanded Stent Length (mm)	Total Drug content on Stent System (µg) (±10%)
EAN22508	2.25	8	42
EAN25008	2.50	8	42
EAN27508	2.75	8	54
EAN30008	3.00	8	54
EAN35008	3.50	8	54
EAN40008	4.00	8	66
EAN45008	4.50	8	66
EAN50008	5.00	8	66
EAN22512	2.25	12	63
EAN25012	2.50	12	63
EAN27512	2.75	12	81
EAN30012	3.00	12	81
EAN35012	3.50	12	81
EAN40012	4.00	12	99
EAN45012	4.50	12	99
EAN50012	5.00	12	99
EAN22516	2.25	16	84
EAN25016	2.50	16	84
EAN27516	2.75	16	108
EAN30016	3.00	16	108

EAN35016	3.50	16	108
EAN40016	4.00	16	132
EAN45016	4.50	16	132
EAN50016	5.00	16	132
EAN22520	2.25	20	105
EAN25020	2.50	20	105
EAN27520	2.75	20	135
EAN30020	3.00	20	135
EAN35020	3.50	20	135
EAN40020	4.00	20	165
EAN45020	4.50	20	165
EAN50020	5.00	20	165
EAN22524	2.25	24	126
EAN25024	2.50	24	126
EAN27524	2.75	24	162
EAN30024	3.00	24	162
EAN35024	3.50	24	162
EAN40024	4.00	24	198
EAN22528	2.25	28	147
EAN25028	2.50	28	147
EAN27528	2.75	28	189
EAN30028	3.00	28	189
EAN35028	3.50	28	189
EAN40028	4.00	28	231
EAN22532	2.25	32	168
EAN25032	2.50	32	168
EAN27532	2.75	32	216
EAN30032	3.00	32	216
EAN35032	3.50	32	216
EAN40032	4.00	32	264
EAN22536	2.25	36	189
EAN25036	2.50	36	189
EAN27536	2.75	36	243
EAN30036	3.00	36	243
EAN35036	3.50	36	243
EAN40036	4.00	36	297
EAN22540	2.25	40	210
EAN25040	2.50	40	210
EAN27540	2.75	40	269
EAN30040	3.00	40	269
EAN35040	3.50	40	269
EAN40040	4.00	40	330
EAN25044	2.50	44	236
EAN30044	3.00	44	306
EAN35044	3.50	44	306
EAN40044	4.00	44	380
EAN25048	2.50	48	256
EAN30048	3.00	48	342
EAN35048	3.50	48	342
EAN40048	4.00	48	414
EAN25052	2.50	52	276
EAN30052	3.00	52	378
EAN35052	3.50	52	378
EAN40052	4.00	52	446

## 2.0 INDICATIONS

The AbLuminus DES+ Sirolimus Eluting Coronary Stent System is indicated for improving luminal diameter in patients with diabetes mellitus, acute myocardial infarction, symptomatic heart disease, stable angina, unstable angina or documented silent ischemia due to atherosclerotic lesions in native coronary arteries with reference vessel diameter 2.25mm to 5.00mm and lesions of length  $\leq$  48 mm.

## 3.0 CONTRAINDICATIONS

The AbLuminus DES+ Stent is contraindicated for use in:

### Condition related to patient:

- Patients with known hyper sensitivity to (i) Sirolimus drug or its structurally related compounds and/or above mentioned biodegradable polymers, (ii) L605 cobalt chromium alloy and (iii) contrast agents that cannot be controlled prophylactically prior to AbLuminus DES+ Stent implantation.

- Arterial spasm
- Patients with an ejection fraction < 30%
- Patients having experienced cardiogenic shock.
- Patients with bleeding diathesis or other disorder e.g. Peptic ulceration or recent cerebrovascular accident, limiting the use of antiplatelet and/or anticoagulation therapy.
- Off-label use (i.e. outside of the approved indications for use). Patient out comes may not be the same as the results observed in clinical studies.

### Condition related to lesion:

- Reference vessel diameter at the lesion site either <2.25 mm or >5.0 mm.
- Heavily calcified lesion not treated with plaque modification strategy.
- Stenosis which cannot be pre-dilated with an angioplasty balloon to a mean luminal diameter of 2.25 mm.

## 4.0 ANTIPLATELET REGIMEN

Continuation of combination treatment with a P2Y12 inhibitor and aspirin after PCI appears to reduce major adverse cardiac events. Patients should receive P2Y12 inhibitor in addition to aspirin for a minimum of 3 months (stable coronary artery disease) or 6 months (acute coronary syndromes) unless there are high bleeding risks. Aspirin (81 mg daily) should usually be continued indefinitely in patients with coronary artery disease.

The optimal duration of dual antiplatelet therapy (DAPT) following DES implantation is unknown, and DES thrombosis may still occur despite continued therapy. Consistent with 2016 ACC/AHA guidelines, longer duration of DAPT may be considered in patients at higher ischemic risk with lower bleeding risk. Decisions about the duration of DAPT are best made on an individual basis and should integrate clinical judgment, assessment of the benefit/risk ratio, and patient preference as per choice of physician. Premature discontinuation or interruption of prescribed antiplatelet medication could result in a higher risk of stent thrombosis, MI or death.

## 5.0 WARNINGS

- Ensure that the inner Tyvek package has not been opened or damaged as this may indicate the sterile barrier has been breached.
- The use of product carries the risks associated with coronary artery stenting, including stent thrombosis, vascular complications, and/or bleeding events.
- Patients with contraindicated/resistant antiplatelet or anticoagulant therapy is a clear call of the treating physician and stenting to be done with relative risk calculation to the patient by the physician.
- Only cardiologists who have received proper training shall perform the stent implantation.
- Ensure that the medical team is trained on the products and their reference system to avoid any error in choosing equipment.
- PTCA in patients who are not acceptable candidates for coronary artery bypass graft surgery requires careful consideration including possible hemodynamic support during the procedure, as treatment of this patient population came special risk.
- In patients with lesion(s) that prevents complete inflation of an angioplasty balloon or proper placement of stent or delivery catheter, stenting should be done post proper vessel bed preparation using various tools like cutting balloon/scoring balloon, rotablation catheter, lithoplasty balloon etc.
- Subsequent restenosis may require repeat dilatation of the arterial segment containing the stent. Long-term outcome following repeat dilatation of endothelialised stents is not known at present.
- This stent deployment device must not be reused in another procedure.
- The performance characteristics of the balloon are degraded during use.
- Direct Stenting is not recommended.
- Implanting a stent may lead to a flow limiting / non-flow limiting dissection in the vessel distal and / or proximal to the stent and may cause acute closure of the vessel requiring additional intervention(CABG, further dilatation, placement of additional stents, or other). Do not expand the stent if it is not properly positioned in the vessel. (See section: Precautions Stent System Removal Precautions see section 6.4.)
- This product is not intended or approved for use in peripheral

applications.

#### **Carcinogenicity, Genotoxicity and Reproductive Toxicology**

The carcinogenicity, and reproductive toxicity of Sirolimus Eluting Coronary Stent System have not been evaluated. However, carcinogenicity, and reproductive toxicity of sirolimus have been investigated in bacterial and mammalian cells in vitro and in laboratory animals in vivo.

#### **Carcinogenicity**

Carcinogenicity studies were conducted in mice and rats. In an 86-week female mouse study at sirolimus doses 30 to 120 times higher than the 2 mg daily clinical dose (adjusted for body surface area), there was a statistically significant increase in malignant lymphoma at all dose levels compared with controls. In a second mouse study at dosages that were approximately 3 to 16 times the clinical dose (adjusted for body surface area), hepatocellular adenoma and carcinoma in males were considered sirolimus-related. In the 104-week rat study at dosages equal to or lower than the clinical dose of 2 mg daily (adjusted for body surface area), there were no significant findings.

#### **Reproductive Toxicity**

The study of the effects of sirolimus on the different stages of the reproductive process were conducted in rats. Fertility was diminished slightly in both male and female rats following oral administration of sirolimus at doses approximately 10 times or 2 times, respectively, the clinical dose of 2 mg daily (adjusted for body surface area). Male rats showed decreased fertility and testicular atrophy (giant cells and hypospermia in the testes and epididymides). In female rats, reduced size of ovaries and uteri was observed.

#### **Genotoxicity**

Genotoxicity test has been performed as per ISO 10993-3:2014 for sirolimus Eluting stent system. The results show that the product doesn't have genetic toxicity in CHO cell line and it is non-mutagenic in salmonella typhimurium Strains TA 1535, TA97a, TA98, TA100 & TA102.

### **6.0 PRECAUTIONS**

#### **6.1 DRUG INTERACTIONS**

The Drug interaction studies have not been conducted with Ablumimus DES+ Stent. Sirolimus is extensively metabolized by cytochrome P4503A4 in the gut wall and liver and under goes efflux from enterocytes of the small intestine by P-glycoprotein(P-gp). Therefore, absorption and subsequent elimination of systemically absorbed drug may be in fluenced by other drugs that affect these proteins. Drugs that may increase Sirolimus blood concentrations include: Calcium channel blockers (nicardipine), Antifungal agents (clotrimazole, fluconazole, itraconazole), Macrolideantibiotics-(clarithromycin, troleandomycin), Gastrointestinal prokinetic agents (cisapride, metoclopramide), other drugs (bromocriptine, cimetidine, danazol, HIV-protease inhibitors). Drugs that may decrease Sirolimus levels include Anticonvulsants (carbamazepine, phenobarbital, phenytoin) and Antibiotics (rifabutin, rifampentine). These lists are not all inclusive and care should be exercised when drugs or other substances that are metabolized by Cytochrome P4503A4 (CYP3A4) are administered concomitantly with implantation of Ablumimus DES+ Stent. The mechanism or mechanisms by which Ablumimus DES+ Stent affects neointimal proliferation has not been completely established. It is however known, that Sirolimus inhibits T-lymphocyte activation and smooth muscle and endothelial cell proliferation in response to cytokine and growth factor stimulation. In cells, it binds to the immunophilin, cytosolic protein FK-binding protein 12 (FKBP12). The Sirolimus - FKBP - 12 complex binds to and inhibits the activation of the mammalian Target of Rapamycin (mTOR), leading to inhibition of cell cycle progression from G<sub>1</sub> to S phase.

#### **6.2 STENT HANDLING PRECAUTIONS**

- For single use only. Do not re-sterilize, reuse or use product reached or exceeded its "Use By" as indicated on product label.
- Do not use if packaging is damaged or opened. The sterility and stability of the product cannot be guaranteed once the Tyvek pouch has been opened. Un-used products should be discarded or returned to Concept Medical B.V. and should not be re-stocked.
- Do not use if stent coating is subjected to abrasions beyond those of normal insertion and delivery.

- Do not rub or scrape the stent coating. Do not use if stent is exposed to abnormal rubbing or contact with objects other than the guide catheter or opened hemostasis valve prior to implantation.
- Do not remove the stent from the delivery balloon as removal may damage the stent and coating and/or lead to stent dislodgement. The Ablumimus DES+ Stent is intended to perform as a whole system. Prior to angioplasty, carefully examine all equipment's / ancillary devices to be used during the procedure.
- The delivery system should not be used in conjunction with other stents.
- Special care must be taken not to touch or in any way disrupt the stent on the balloon. This is most important during stent system removal from packaging, placement over guidewire, and advancement through rotating hemostatic valve adaptor and guiding catheter hub.
- Do not "roll" the mounted stent with your fingers as this action may loosen the stent from the delivery balloon and can cause subsequent dislodgment, or loss of drug coating.
- Use only the appropriate balloon inflation media. Do not use air or any gaseous medium to inflate the balloon as this may cause uneven expansion and difficulty in deployment of the stent.
- Do not attempt to straighten the proximal shaft (hypotube) as it may cause the catheter to break if it is accidentally bent.
- In the event the Ablumimus DES+ Stent is not successfully deployed, follow product return procedures and avoid handling of the stent with bare hands.

### **6.3 STENT PLACEMENT PRECAUTIONS**

- **Do not prepare, introduce negative pressure, or pre-inflate the delivery system prior to stent deployment other than as directed.**
- All supporting accessories used during use of Ablumimus DES+ Sirolimus eluting coronary stent system shall be undamaged and used before expiry.
- Use balloon purging technique described in section: Delivery System Preparation see section 9.3.
- The labelled stent diameter refers to expanded stent inner diameter.
- **The vessel should be pre-dilated with an appropriately sized balloon.**
- When treating multiple lesions, the distal lesion should be initially stented, followed by stenting of the proximal lesion. Stenting in this order obviates the need to cross the proximal stent in placement of the distal stent and reduces the chances for dislodging the proximal stent.
- **Guiding Catheter and Guidewire selection**
- Only guiding catheters and guidewires indicated for use in coronary angioplasty should be used. Use a guiding catheter with a minimum internal diameter of 0.058" (1.47mm).
- Guidewire diameter recommended is 0.014" (0.35 mm)
- **Stent Selection**
- Careful stent sizing is important for successful stenting. In general, the stent size should match the diameter of the reference vessel and the length of the lesion to be stented. The diameter of the expanded delivery system must not be superior to the segment diameter, proximal and distal to the lesions. Length of stent is selected 2 mm at least longer distally and proximally from the lesion length or affected segment.
- **Use of Multiple Stents:** When multiple drug-eluting stents are required to cover the long lesion, use additional stents as necessary to adequately cover the lesion. The Sirolimus eluting coronary stent system clinical program data recommends that same stent shall be used for overlapping in terms of drug and stent alloy. Use of another stent shall include but not limited to identical in chemical composition of overlapping stent physicochemical properties of the surface (especially with reference to surface treatment of the metals).
- The possible interactions of Ablumimus DES+ Stent with other drug-eluting or coated stents have not been evaluated and should be avoided whenever possible.
- When using two wires, care should be taken when introducing, torquing and removing one or both guide wires to avoid entanglement. In this situation, it is recommended that one guide wire be completely withdrawn from the patient before removing any additional equipment.
- Do not expand or proceed to deploy the stent if it is not properly positioned in the vessel. (See section 6.4: Stent System Removal Precautions).

- Placement of a stent has the potential to compromise side branch patency. Balloon pressures should be monitored during inflation. **Do not exceed rated burst pressure as indicated on the product label.** (Refer to labelling) Use of pressures higher than specified on the product label may result in a ruptured balloon with possible intimal damage and dissection.
- **Do not attempt to pull an unexpanded stent back through the guiding catheter, as dislodgement of the stent from the balloon may occur.** Remove as a single unit. (See section 6.4: Stent System Removal Precautions).
- Stent retrieval methods (use of additional wires, snares and/or forceps) may result in additional trauma to the coronary vasculature and/or the vascular access site. Complications may include bleeding, hematoma or Pseudo aneurysm. Administration of appropriate anticoagulant, antiplatelet and coronary vasodilator therapy is critical to successful stent implantation.

#### 6.4 STENT SYSTEM REMOVAL PRECAUTIONS

Should unusual resistance be felt at any time during either lesion access or removal of the Delivery System prior to stent implantation, the entire system should be removed as a single unit.

##### When removing the Delivery System as a single unit:

- Do not attempt to pull an unexpanded stent back through the guiding catheter while engaged in the coronary arteries, as stent damage or stent dislodgement from the balloon may occur. Advance the guidewire into the coronary anatomy as far distally as safely possible.
- Position the proximal balloon marker just distal to the tip of the guiding catheter.  
Note: If it is necessary to maintain guidewire position, then it must either be converted to an exchange wire length or a second guidewire must be inserted.
- Tighten the rotating hemostatic valve to secure the Delivery System to the guiding catheter. Remove the guiding catheter and Delivery System as a single unit.

Note: After use, this product may be a biohazard. Handle and dispose of all such devices in accordance with accepted medical practice and applicable government regulations.

Do not attempt to pull the guiding catheter and delivery catheter through the femoral sheath. When the distal tip of the guiding catheter reaches the distal end of the femoral sheath, remove sheath, guiding catheter, and delivery system as a single unit and replace sheath per hospital protocol.

Failure to follow these steps and / or applying excessive force to the Delivery System can potentially result in stent dislodgement or damage to the stent and / or Delivery System components.

If it is necessary to retain guidewire position for subsequent artery / lesion access, leave the guide wire in place and remove all other system components.

#### 6.5 POST PROCEDURE PRECAUTIONS

Great care must be exercised when crossing a newly deployed stent with a coronary guidewire, IVUS catheter, balloon or another stent delivery system to avoid disrupting the stent geometry.

#### 6.6 MAGNETIC RESONANCE (MR) IMAGING

Non-clinical testing has demonstrated that Sirolimus Eluting Coronary Stent System in single and in overlapped configurations up to 70 mm in length is MR Conditional. It can be safely scanned under the following conditions:

- Static magnetic field of 1.5-Tesla (1.5 T) or 3-Tesla (3 T)
  - Maximum spatial field gradient of 2,630 Gauss/cm (26.30 T/m)
  - Maximum whole body averaged specific absorption rate (SAR) of:
    - 2.0 W/kg (Normal Operating Mode) at 1.5 T
    - 2.0 W/kg (Normal Operating Mode) at 3 T
- Sirolimus Eluting Coronary Stent System should not migrate in this MRI environment. Non clinical testing at field strengths greater than 3 Tesla has not been performed to evaluate stent migration or heating. MRI at 1.5 or 3 Tesla may be performed immediately following the implantation of Sirolimus Eluting Coronary Stent System.

**RF heating:** Sirolimus Eluting Coronary Stent System single and in overlapped configurations up to 70 mm in length, are expected to produce a maximum temperature rise of less than 6°C after 15 minutes of continuous scanning at an MR system reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg.

**MR Artifact:** In non-clinical testing, the image artifact caused by the device extends approximately 0.8 cm from Sirolimus Eluting Coronary Stent System when imaged with gradient-echo and spin echo pulse sequences in a 3 T MRI system.

#### 7.0 INDIVIDUALIZATION OF TREATMENT

See also section 8.0: Use in Special Populations.

The risks and benefits of drug eluting stents should be considered for each patient before use of the Ablumimus DES+ Stent. Cardiologists are responsible for assessing patient appropriateness for stent implantation prior to procedure.

#### 8.0 USES IN SPECIAL POPULATION

**PREGNANCY:** There are no data available for use of Ablumimus DES+ stent in pregnant women.

**LACTATION:** The effects of Sirolimus drug during lactation have not been evaluated.

**GENDER:** Clinical studies of Sirolimus based stents did not find any significant differences in safety and effectiveness for male and female patients.

**ETHNICITY:** Clinical studies of the Ablumimus DES+ stent did not include sufficient numbers of patients to assess for differences in safety and effectiveness due to ethnicity.

**PEDIATRIC USE:** The safety and effectiveness of Ablumimus DES+ Stent below the age of 18 years have not been established.

Carefully consider whether it is appropriate to use this stent in the above patient populations.

The safety and effectiveness beyond five years or of using mechanism atherectomy devices (directional atherectomy catheters, rotation atherectomy catheters) or laser angioplasty catheters to treat in-stent stenosis has not been established.

#### 9.0 OPERATOR'S MANUAL

##### 9.1 INSPECTION PRIOR TO USE

- Inspect the stent delivery system package for damage to the sterile barrier.
- Carefully remove the system from the package and inspect the delivery catheter for bends, kinks and other damage.
- Remove stylet (protection wire) carefully.
- Inspect the stent to ensure that it has not been damaged or displaced from its original position on the balloon. Verify that the stent is positioned between the proximal and distal balloon markers.
- Note the position of the stent relative to the proximal and distal marker bands for use as reference under fluoroscopy.
- Do not use device if any defects are noted.

##### 9.2 MATERIALS REQUIRED

Quantity	Material
1	Appropriate guiding catheter(s) with minimum inner diameter of 5F (~ 0.058")
1	Pre-dilatation balloon catheter
1	10-20cc syringe
1000IU	Heparin per 500cc Normal Saline (HepNS)
1	0.014" guidewire ≥ 175 cm
1	Rotating hemostatic valve
N/A	Contrast diluted 1:1 with normal saline
1	Inflation device
1	3-way Stopcock
1	Torque device

##### 9.3 DELIVERY SYSTEM PREPARATION

1. Prepare the inflation device or syringe with diluted contrast medium.
2. Attach the inflation device or syringe to the 3-way stopcock; attach to balloon inflation port hub.  
Caution: Do not apply negative or positive pressure to balloon at this time as this can cause premature dislodgement of the stent.
3. Open the stop cock to stent delivery system.
4. Leave the inflation device or syringe on neutral.

##### 9.4 STENT DELIVERY PROCEDURE

1. Prepare the vascular access site according to standard PTCA practice.
2. Pre-dilate the lesion with a balloon diameter 0.5 mm smaller than the stent and a balloon length equal to or shorter than the target lesion length and shorter than the length of the stent to be implanted. Immediately prior to back loading the stent

delivery catheter on to the guidewire, flush the guidewire lumen with HepNS according to hospital protocol. Avoid contact with stent.

**NOTE:** Stenting should always be done followed by proper vessel bed preparation. Stent contact with fluid has the possibility of initiating drug release. Fluid contact time should be limited to immediately prior to loading the delivery catheter on the guidewire. When flushing the delivery catheter, ensure that damage to the stent does not occur.

3. Back load the stent delivery system onto the proximal portion of the guide wire while maintaining its position across the target lesion.
4. Open rotating / Push-Pull hemostatic valve on the guiding catheter hub as widely as possible and close when the stent has been advanced.
5. Advance the stent delivery system over the guidewire to the target lesion under fluoroscopic guidance. Use the radiopaque balloon markers to position the stent across the lesion. Perform angiography to confirm the position of the stent.  
**NOTE:** If resistance is felt, do not force passage. Resistance may indicate a problem and may result in damage to the vessel, stent or stent dislodgement if it is forced. Remove the stent delivery system and guiding catheter as a single unit. (See section 6.4: Stent System Removal Precautions)

#### 9.5 DEPLOYMENT PROCEDURE

1. Before deployment, reconfirm the correct position of the stent relative to the target lesion via the radiopaque balloon markers.
2. Ensure that the 3-way stopcock on the stent delivery system is open to the inflation device and apply negative pressure to purge the balloon air.
3. Turn the 3-way stopcock on the stent delivery catheter off to the balloon port and purge the inflation device of air. Open the side port of the 3-way stopcock to the delivery system.
4. Under fluoroscopic visualization, inflate the balloon to at least the nominal pressure to deploy the stent, but do not exceed the labelled rated burst pressure (RBP). It is highly recommended that for optimal deployment of stent and release of the drug from exposed part of balloon, it is advised to keep the stent delivery system inflated for 45 seconds. However, it is at the sole discretion of the doctor and also on patient tolerance. In case needed, multiple inflations can be done for optimal deployment of the stent and drug delivery from exposed parts of balloon. Optimal expansion requires the stent to be in full contact with the artery wall, with the stent internal diameter matching the size of the reference vessel diameter. Ensure that the stent is not under dilated.
5. Deflate the balloon by pulling a vacuum with the inflation device for upto 15 seconds. Make sure the balloon is fully deflated before attempting any movement of the system.
6. Confirm adequate stent expansion and balloon deflation by angiographic injection through the guiding catheter.
7. If more than one Abluminus DES+ Stent is needed to cover the lesion and balloon treated area, adequately overlap the stents.

#### 9.6 REMOVAL PROCEDURE

1. Ensure that the balloon is fully deflated.
2. Fully open the rotating or push pull hemostatic valve.
3. While maintaining the guidewire position and negative pressure on the inflation device, withdraw the stent delivery system.
4. Tighten rotating or push pull hemostatic valve.
5. Should unusual resistance be felt at anytime during removal of stent delivery system after stent implantation, the entire system should be removed as a single unit.
6. Repeat angiography to assess the stented area. If adequate expansion has not been obtained, exchange back to original stent delivery catheter or exchange to another balloon of appropriate balloon diameter to achieve proper stent apposition to the vessel wall.

#### 9.7 FURTHER DILATATION OF STENTED SEGMENTS

1. If an adequate expansion has not been obtained, re-advance the stent delivery system or exchange for another balloon catheter of appropriate balloon diameter to achieve proper stent apposition to the vessel wall.  
**NOTE:** Post-dilatation should be performed with in the stented segment. Do not dilate beyond the stent edges.
2. Reconfirm stent position and post procedure angiographic result. Repeat the inflations until optimal stent deployment and TIMI-III flow is achieved. Final stent diameter should match reference vessel diameter.

## 10.0 POTENTIAL ADVERSE EVENTS

### POTENTIAL ADVERSE EVENTS RELATED TO SIROLIMUS

The side effects/complications may be associated with, but not limited to, the use of Sirolimus:

- Anemia, Diarrhea, Dry Skin, Headache, Hematuria, Infection, Injection site reaction, Pain (abdominal, arthralgia, injection site), Rash

### POTENTIAL ADVERSE EVENTS RELATED TO POLYMER

The side effects/ complications of the polymer are not different than those of other commercially available drug eluting stent coatings and may include, but are not limited to the following:

- Allergic reaction
- Focal inflammation at the site of stent implantation
- Restenosis of the stented artery

### POTENTIAL RISKS ASSOCIATED WITH PERCUTANEOUS CORONARY DIAGNOSTIC AND TREATMENT PROCEDURES

Other risks associated with using this device are those associated with percutaneous coronary diagnostic (including angiography and IVUS) and treatment procedures. These risks (in alphabetical order) may include but are not limited to the following:

- Abrupt vessel closure
- Allergic reaction (to contrast, antiplatelet therapy, stent material, or drug and polymer coating)
- Arrhythmias, including ventricular fibrillation and ventricular tachycardia
- Balloon rupture
- Bleeding
- Cardiac tamponade
- Coronary artery occlusion, perforation, rupture, or dissection
- Coronary artery spasm
- Coronary Aneurysm
- Death
- Embolism (air, tissue, device, or thrombus)
- Emergency surgery: peripheral vascular or coronary bypass.
- Failure to deliver the stent
- Hematoma
- Hemorrhage requiring transfusion
- Hypotension / hyper tension
- Incomplete stent apposition
- Myocardial infarction (MI)
- Pyrexia
- Pericarditis
- Peripheral ischemia / peripheral nerve injury
- Renal Failure
- Respiratory Failure
- Restenosis of the stented artery
- Shock / pulmonary edema
- Stable or Unstable angina
- Stent deformation, collapse, or fracture
- Stent migration or embolization
- Stent misplacement
- Stroke / transient ischemic attack
- Thrombosis (acute, sub acute or late)

## 11.0 HOW SUPPLIED

**Sterile:** This device is sterilized with Ethylene Oxide Gas and is non-pyrogenic.

**Contents:** One Abluminus DES+ sirolimus eluting stent mounted on Rapid exchange delivery system.

**Storage:** Store between 8°C to 25°C at dry place and away from direct sunlight.

**Disposal:** Dispose of device shall be in accordance with guidelines of local regulations.

## 12.0 COMPLIANCE CHART, INFLATION PRESSURE

Refer to labeling.

## 13.0 IMPLANT CARD

In addition to these Instructions for Use booklet, a Product Implant Card is provided in the package. It is recommended that patients keep this card in their possession for procedure/stent identification.

## 14.0 DISCLAIMER OF WARRANTY

Concept Medical B.V. does not provide any expressed or implied warranty, including without limitation, any implied warranty of merchantability, fitness for a particular purpose or non-infringement, with respect to the product(s) described in this publication and Concept Medical B.V. hereby expressly disclaims any such warranties.

Descriptions or specifications in Concept Medical B.V. printed matter, including this publication, are meant solely to generally describe the relevant product at the time of manufacture and do not constitute any express warranties.

No person has the authority to bind Concept Medical B.V. to any representation or warranty.

Under no circumstances shall Concept Medical B.V. be liable for any medical or other expenses or direct, indirect, incidental, consequential, punitive or exemplary damages with respect to the products described herein, including without limitation, in respect of any use, defect, failure or malfunction, whether based upon warranty, contract, tort or other legal theory, even if Concept Medical B.V. has been advised or should have been aware of the possibility of such damages, other than as expressly provided by specific law without possibility of exclusion.

The exclusions, disclaimers and limitations contained in this section shall be interpreted in such manner as to be effective and valid under applicable law, but in case any one or more of the provisions contained herein shall, for any reason, be held to be invalid, illegal or unenforceable in any respect, such provision shall be ineffective to the extent, but only to the extent, of such invalidity, illegality or unenforceability without invalidating the remainder of this section or any other provisions hereof.

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Nijedna osoba nema ovlašćenje da zavezuje Concept Medical B.V. za bilo koje predstavljanje ili garanciju.