 **KONAR-MF™**

Instructions for Use	
KONAR-MF™ VSD Occluder Instructions for Use	EN 1-6
شاميلعت	
مطبق إغلاق ثقب الحاجز البطني KONAR-MF™ تعليمات الاستخدام	AR 7-11
Brugsanvisning	
Brugsanvisning til KONAR-MF™ VSD-okklusjonsanordning	DA 12-17
Gebrauchsanweisung	
KONAR-MF™ VSD Okkluder Gebrauchsanweisung	DE 18-24
Οδηγίες Χρήσης	
Συσκευή Απόφραξης VSD KONAR-MF™ Οδηγίες Χρήσης	EL 25-31
Instrucciones de uso	
Instrucciones de uso del oclusor de VSD KONAR-MF™	ES 32-37
Mode d'emploi	
Système de fermeture d'ASV KONAR-MF™ Mode d'emploi	FR 38-43
Használati utasítás	
KONAR-MF™ VSD-okkluder használati utasítás	HU 44-49
Istruzioni per l'uso	
Occlusore per VSD KONAR-MF™ Istruzioni per l'uso	IT 50-55
Bruksanvisning	
Bruksanvisning for KONAR-MF™ VSD-okklusjonsanordning	NO 56-61
Instrukcja użytkowania	
Okluder KONAR-MF™ VSD Instrukcja obsługi	PL 62-67
Instruções de utilização	
Oclusor KONAR-MF™ VSD Instruções de utilização	.PT 68-73
Инструкция по эксплуатации	
Окклюдер KONAR-MF™ VSD Инструкция по применению	.RU 74-80
Brugsanvisning	
Brugsanvisning för KONAR-MF™ VSD ocklusionsanordning	SV 81-86
Kullanım Talimatları	
KONAR-MF™ VSD Oklüder Kullanım Talimatları	TR 87-92

	<p>Date of manufacture تاريخ التصنيع Fremstillingsdato Herstellungsdatum Ημερομηνία κατασκευής Fecha de fabricación Date de fabrication Gyártás dátuma Data di fabbricazione Produksjonsdato Data produkcji Data de fabrico Дата изготовления Tillverkningsdatum Üretim tarihi</p>
	<p>Manufacturer الشركة المصنعة Producent Hersteller Κατασκευαστής Fabricante Fabricant Gyártó Produttore Produsent Producent Fabricante Производитель Tillverkare Üretici</p>
	<p>Use-by date يُستخدم من قبل Anvendes inden Verwendbar bis Χρήση έως Fecha de caducidad Date de péremption Felhasználhatósági idő Utilizzare entro Utløpsdato Użyć przed Utilizar por Использовать до Förbrukningsdatum Son kullanma tarihi</p>

<div data-bbox="146 220 291 303" style="border: 1px solid black; padding: 5px; text-align: center; width: fit-content; margin: auto;"> SN </div>	<p>Serial number الرقم المسلسل Seriennummer Seriennummer Αριθμός σειράς Número de serie Numéro de série Sorozatszám Numero di serie Seriennummer Numer seryjny Número de série Серийный номер Seriennummer Seri numarası</p>
<div data-bbox="151 614 285 703" style="border: 1px solid black; padding: 5px; text-align: center; width: fit-content; margin: auto;"> LOT </div>	<p>Batch code كود التشغيل Batchkode Chargennummer Κωδικός παρτίδας Código de lote Numéros de lot Sarzs-kód Codice del lotto Batchnummer Kod partii Código do lote Код серии Satskod Parti kodu</p>
<div data-bbox="107 1016 329 1089" style="border: 1px solid black; padding: 5px; text-align: center; width: fit-content; margin: auto;"> STERILE EO </div>	<p>Sterilized using ethylene oxide معقم باستخدام أكسيد الإيثيلين Steriliseret med etylenoxid Mit Ethylenoxid sterilisiert Αποστειρωμένο με αιθυλενοξειδίο Esterilizado con óxido de etileno Stérilisé à l'oxyde d'éthylène Etilénokziddal sterilizálva Sterilizzato con ossido di etilene Steriliseret med etylenoksid Wyjałowiono tlenkiem etylenu Esterilizado utilizando óxido de etileno Стерилизовано этиленоксидом Steriliserad med etylen Etilen oksitle sterilize edilmiştir</p>

	<p>Non-pyrogenic غير مُؤَلِّدٌ لِلْحُمَّى Ikke-pyrogen Nicht pyrogen Μη πυρετογόνο Apirógeno Non pyrogène Nem pirogén Apirogeno Ikke-pyrogent Produkt niepirogenny Não pirogénico Апирогенно Icke-pyrogen Pirojenik deǵildir</p>
	<p>Do not reuse لا يُسْتخدَم ثَانِيَةً Må ikke genanvendes Nicht wiederverwenden Να μην επαναχρησιμοποιείται No reutilizar Ne pas réutiliser Nem újrafelhasználható Non riutilizzare Ikke til gjenbruk Nie używać ponownie Não reutilizar Не использовать повторно Får ej återanvändas Yeniden kullanmayın</p>
	<p>Consult Instructions For Use ارجع إلى تعليمات الاستخدام Se bruksanvisningen Gebrauchsanweisung beachten Συμβουλευτείτε τις Οδηγίες χρήσης Consultar las instrucciones de uso Consulter le mode d'emploi Olvassa el a használati utasítást Consultare le istruzioni per l'uso Se bruksanvisningen Patrz Instrukcja stosowania Consulte as Instruções de utilização Обратитесь к инструкции по применению Se bruksanvisning före användning Kullanım Talimatlarına Başvurun</p>

Keep away from sunlight
يُحفظ بعيدًا عن أشعة الشمس
Holdes væk fra sollys
Von Sonnenlicht fernhalten
Φυλάξτε μακριά από το ηλιακό φως
Mantener alejado de la luz solar
Conserver à l'abri de la lumière
Napfénytől távol tartandó
Tenere lontano dalla luce del sole
Hold unna sollys
Chronić przed promieniowaniem słonecznym
Mantenha afastado da luz solar
Защищать от солнечного света
Får ej utsättas för solljus
Güneş ışığından koruyun

Keep dry
يُحفظ جافًا
Holdes tør
Trocken halten
Να διατηρείται στεγνό
Mantener seco
Conserver au sec
Százazon tartandó
Conservare in luogo asciutto
Hold tørt
Przechowywać w suchym miejscu
Manter seco
Хранить в сухом месте
Förvaras torrt
Kuru ortamda saklayın

Catalogue number
رقم الكتالوج
Katalognummer
Katalognummer
Αριθμός καταλόγου
Número de catálogo
Référéncé
Katalógusszám
Numero di catalogo
Katalognummer
Numer katalogowy
Número do catálogo
Νομep по каталогу
Artikelnummer
Katalog numarası

	<p>Do not use if package is damaged لا تستخدم العبوة في حالة تلفها Må ikke bruges, hvis pakningen er beskadiget Bei beschädigter Verpackung nicht verwenden Να μην χρησιμοποιείται αν υπάρχει ζημιά στη συσκευασία No utilizar si el paquete está dañado Ne pas utiliser si l'emballage est endommagé Sérült csomagolás esetén ne használja Non utilizzare se la confezione è danneggiata Ikke bruk hvis emballasjen er skadet Nie używać, jeśli opakowanie jest uszkodzone Não utilizar se a embalagem estiver danificada Не использовать в случае повреждения упаковки Använd inte om förpackningen är skadad Paket hasar görmüşse kullanmayın</p>
	<p>Do not resterilize لا تعيد تعقيمه Må ikke gensteriliseres Nicht erneut sterilisieren Να μην επαναποστειρώνεται No reesterilizar Ne pas restériliser Ne sterilizálja újra Non risterrilizzare Ikke resteriliser Nie wyjalawiać Não voltar a esterilizar Не стерилизовать повторно Får ej omsteriliseras Tekrar sterilize etmeyin</p>
	<p>Temperature limit حدود درجة الحرارة Temperaturgrænse Temperaturgrenze Όρια θερμοκρασίας Limite de temperatura Limite de température Hőmérséklet-határérték Limite di temperatura Temperaturgrense Dopuszczalna temperatura Limite de temperatura Ограничение температуры Temperaturgräns Sıcaklık limiti</p>

Authorised representative in the European Community
الممثل المعتمد في الاتحاد الأوروبي
Autoriseret repræsentant i Det Europæiske Fællesskab
Bevollmächtigter in der EU
Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Ένωση
Representante autorizado en la comunidad Europea
Représentant agréé dans la communauté européenne
Meghatalmazott képviselő az Európai Közösségben
Rappresentante autorizzato UE
Autorisert representant i EU
Autoryzowany przedstawiciel w UE
Representante autorizado da UE
Уполномоченный представитель в ЕС
Godkänd representant i Europeiska Gemenskapen
AB yetkili temsilcisi

MR conditional
أمان نينالر المغناطيسي طولالمشر
Betinget MR-sikker
Bedingte magnetische Resonanzsicherheit
Ασφαλές για μαγνητική τομογραφία υπό προϋποθέσεις (MR Conditional)
RM condicional
Conditions d'RM
MR-kondicionális
Compatibilità RM
MR-betinget
Warunkowe bezpieczeństwo w badaniach MR
Compatibilidade MR
МРТвыполнятьс соблюдениемограничений
MR-villkor
MR koşullu

	<p>Medical device طبي جهاز Medicinsk udstyr Medizinisches Gerät Ιατροτεχνολογικό προϊόν Producto sanitario Dispositif médical Orvostechnikai eszköz Dispositivo medico Medisinsk utstyr Wyrób medyczny Dispositivo médico Μεδικινσκη οηδελιη Medicinteknisk produkt Tibbi Cihaz</p>
	<p>Importer المستورد Importør Importeur Εισαγωγέας Importador Importateur Importör Importatore Importør Importer Importador Импортер Importör İthalatçı</p>
	<p>Made in china الصين في صنع Fremstillet i Kina Hergestellt in China Κατασκευασμένο στην Κίνα Hecho en China Fabriqué en Chine Kínában készült Fabbricato in Cina Fremstilt i Kina Wyprodukowano w Chinach Fabricado na China Сделано в Китае Tillverkad i Kina Çin'de üretilmiştir</p>

Unique device identifier
الفريد الجهاز مُعرّف
unik udstyrsidentifikation
Eindeutige Geräteerkennung
Αποκλειστική ταυτοποίηση ιατροτεχνολογικού προϊόντος
Identificación única del producto
Identifiant unique du dispositif
Egyedi eszközazonosító
Identificativo univoco del dispositivo
Unik utstyrsidentifikator
Unikalny identyfikator urządzenia
Identificador único do dispositivo (UDI)
Уникальный идентификатор изделия
Unik enhetsidentifikation
Benzersiz cihaz tanımlayıcısı

Double sterile barrier system
مزدوج مُعقم حاجز نظام
System med steril dobbelt barrier
System mit doppelter Sterilbarriere
Σύστημα φραγμού διπλής αποστείρωσης
Sistema de doble barrera estéril
Système de double barrière sterile
Dupla steril gátrendszer
Sistema a doppia barriera sterile
Dobbelt sterilt barrieresystem
System podwójnej bariery sterylnej
Sistema de barreira estéril dupla
Двойная барьерная система для стерилизации
Dubbla sterila barriärsystem
Çifte steril bariyer sistemi

1. Device Description

- The KONAR-MF™ VSD Occluder is a percutaneous, transcatheter closure device for the non-surgical closure of Ventricular Septal Defects (VSDs) in patients. The occluder is a self-expandable, double-disc device made from a Nitinol wire mesh. The two discs are linked together by a cone-shaped waist corresponding to the size of the VSD. The waist of the four largest models are sewn with PTFE membranes securely using Nylon threads in order to increase its occlusion ability and reduce the residual shunts, while the four smallest models have no membrane in it.
- The KONAR-MF™ VSD Occluder can be used in combination with the SteerEase™ Introducer, which consists of a coil reinforced sheath, a dilator, a loader, a haemostatic valve and a delivery cable. The introducer is used to advance the KONAR-MF™ VSD Occluder to the proper position. When the occluder is released from the sheath, one disc expands on each side of the defect, and the expanding waist closes the VSD tunnel in the septum between the left and right ventricles.

2. Indications and Usage

- The KONAR-MF™ VSD Occluder is a percutaneous, transcatheter closure device intended for the occlusion of VSDs.
- Indications:
 - ◇ Patients have VSDs, which effects the cardiac hemodynamics.
 - ◇ Patients older than 6 months and weigh more than 8kg.
 - ◇ Isolated perimembranous VSD with or without Membranous Septal Aneurysm (MSA).
 - ◇ Muscular VSD.
 - ◇ Size of the VSD is larger than or equal to 2mm, less than 10mm.
 - ◇ Upper margin of VSD to aortic valve distance >2 mm for models 6-4, 8-6, 10-8, 12-10, 14-12 and >2.5mm for models 5-3, 7-5, 9-7, in case of PmVSD.
 - ◇ Patient does not present with an aortic valve prolapsing into the VSD and patient does not exhibit aortic valve regurgitation.

3. Contraindications

- Patient has extensive congenital cardinal anomaly which can only be adequately repaired by cardiac surgery.
- Presence of thrombus at the intended site of implant, or documented evidence of venous thrombus in the vessels through which access to the defect is gained.
- Active endocarditis or other infections producing bacteria.
- The implant of KONAR-MF™ VSD Occluder would cause an obvious interference with the aortic valve or the atrioventricular valve.
- Patients with severely increased pulmonary vascular resistance and a right-to-left shunt and patients with documented irreversible pulmonary vascular disease.
- Patients with contraindications to anti-platelet therapy or agents.

4. Warnings

- Patient allergic to nickel may suffer an allergic reaction to the device.
- The KONAR-MF™ VSD Occluder and Introducer should only be used by those physicians trained in transcatheter defect closure techniques.
- Physicians must be prepared to deal with emergency situations, which require removal of embolized devices that result in critical hemodynamic compromise. This includes the availability of an onsite surgeon.

- Embolized devices must be removed. Embolized devices should not be withdrawn through intracardiac structures unless they have been adequately collapsed within a sheath.
- For single use only. Do not reuse or resterilize. Structural integrity and/or function may be impaired or be lost through cleaning, resterilization, or reuse and may cause adverse patient reactions. Lifetech will not be responsible for any direct or consequential damages or expenses resulting from reuse of any of the components in the KONAR-MF™ VSD Occluder.
- Do not use, if the sterile barrier has been compromised in any way.
- Do not release the KONAR-MF™ VSD Occluder from the delivery cable if the device does not conform to its original configuration or if the device position is unstable. Recapture the device and redeploy. If still unsatisfactory, recapture the device and replace with a new device. Implantation of this device may not eliminate or substitute the need for Coumadin in patients with VSD and paradoxical embolism.
- Device closure in patients who have suffered a previous thromboembolic stroke should be discussed with the patient or family. In addition, consultation with a neurologist and hematologist is suggested to determine if the benefit of device closure outweighs the risk.
- Double check the connection safety before loading the device into the loader.
- Retrieve the device when it is hard to release the device from the delivery cable.
- When III° atrioventricular block occurs during VSD closure operation, stop implanting and send for surgical treatment.
- If II° atrioventricular block occurs before the occluder placement, stop to observe for three to five minutes. Implantation can go on only if the atrioventricular block disappears after observation. If II° atrioventricular block occurs after the occluder placement, the patient should be administered with steroid cortex hormone (e.g. Dexamethasone 0.1mg/kg/T). The occluder can be released if the atrioventricular block disappears in fifteen minutes observation, otherwise withdraw the device and send for surgical treatment.

5. Precautions

- The KONAR-MF™ VSD Occluder is single use device. The production date and expiry date are marked on the label. Do not reuse or resterilize them.
- Keep away from sunlight, keep dry, store between 10 °C-30 °C.
- Pre-procedure, procedure and post-procedure medication.
 - ◇ Pre-procedure

Aspirin ([3-5 mg]/kg/day) is to be started at least 24 hours prior to the procedure. Anti-platelet and anticoagulation medication is optional if the patient has contraindication to aspirin.
 - ◇ Procedure

Patient should be fully heparinized throughout the procedure with a minimum active clotting time (ACT) of 250 seconds prior to device insertion.

Transesophageal echocardiography (TEE) or similar imaging equipment is recommended as an aid in placing the Occluder.
- Post-procedure

Patients should take appropriate endocarditic prophylaxis for 6 months following device implantation. The decision to continue endocarditic prophylaxis beyond 6 months is at the discretion of the physician. Patients, who showed any residual shunt, should be inspected by echocardiographic examination for the residual shunt until the defect is closed completely.
- MR conditional

A patient with an implanted KONAR-MF™ VSD Occluder can be scanned safely immediately after

placement of the device under the following conditions:

- ◇ Static magnetic field of 3 T or less.
- ◇ Spatial gradient magnetic field of 720 G/cm or less.
- ◇ Maximum MR system-reported, whole-body-averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning.

Remark: MRI image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the device. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

- Use in Specific population

- ◇ Pregnancy- Care should be taken to minimize the radiation exposure to the fetus and the mother.
- ◇ Nursing mothers- There has been no quantitative assessment of the presence of leachable in breast milk.

6. Potential Adverse Events

- Placement of the KONAR-MF™ VSD Occluder involves using standard interventional cardiac catheterization techniques. The following adverse events might be expected from interventional cardiac catheterization techniques.

- | | |
|--------------------------------------|---|
| ◇ Allergic reaction | ◇ Arrhythmias |
| ◇ Air embolus | ◇ Thrombus formation |
| ◇ Septal tearing | ◇ Device migration |
| ◇ Injury to the nerve or vessel wall | ◇ Cardiac perforation |
| ◇ Incomplete sealing of the defect | ◇ Infection |
| ◇ Hemolysis | ◇ Cardiac tamponade |
| ◇ III° Atrioventricular block | ◇ I° and II° Atrioventricular block |
| ◇ Aortic regurgitation | ◇ Left ventricular enlargement or dysfunction |
| ◇ Fever | ◇ TIA |
| ◇ Residual shunt | ◇ Apnea |
| ◇ Excessive bleeding | ◇ endocarditis |
| ◇ Hematom | ◇ Pericardial effusion |

7. Product features

- The occluders are available in the following (Figure 1 & Table 1).



Figure 1 The KONAR-MF™ VSD Occluder**Table 1** List of specifications of the KONAR-MF™ VSD Occluder and recommended introducers

Unit: mm

Occluder Specification	ΦD	ΦD1	ΦD2	L	Recommended introducer specification
LT-MFO-5-3	10	3	5	4	4-5F
LT-MFO-6-4	10	4	6		4-5F
LT-MFO-7-5	12	5	7		4-5F
LT-MFO-8-6	12	6	8		4-5F
LT-MFO-9-7	14	7	9		6F
LT-MFO-10-8	14	8	10		6F
LT-MFO-12-10	16	10	12		7F
LT-MFO-14-12	18	12	14		7F

8. Directions for Use

- The SteerEase™ Introducer (Refer to Figure 2)
 - ◇ Sheath: The sheath is used to advance the devices to its desired position within the heart.
 - ◇ Hemostatic valve: The haemostatic valve at the proximal end of the sheath minimizes bleeding. The sideport with the flexible extension tube and stopcock is used to flush the system.
 - ◇ Dilator: The dilator is used to ease penetration of the tissues and the vessel wall
 - ◇ Loader: The loader is used to introduce the occluder with the attached delivery cable into the sheath.
 - ◇ Delivery cable: The delivery cable is used to advance the occluder through the sheath, holding it in position while the sheath is pulled back to deploy the occluder. The handle and delivery cable are also used to retrieve and/or reposition the occluder if the size, position or the expansion of the occluder is deemed to be unsatisfactory.
 - ◇ Handle: The plastic vise screwed top the proximal end of the delivery cable is used to facilitate direction control and serve as “release control handle” for disconnecting (releasing) the occluder from the delivery cable.

**Figure 2** The KONAR-MF™ VSD Occluder and the SteerEase™ Introducer

- Directions for Use
 - ◇ Perform procedure under either general or conscious sedation.
 - ◇ Give patient a dose of an appropriate antibiotic during the catheterization procedure.
 - ◇ Obtain access to the femoral artery, femoral vein, and/or right internal jugular vein.
 - ◇ Administer heparin to achieve an activated clotting time of greater than 250 seconds throughout the procedure.
 - ◇ Perform standard right and left heart catheterization. Assess the pulmonary vascular resistance.


Figure 3 Illustration of loading the device into a loader

- ◇ Measure the size of the defect and its distance to the aortic valve by TEE and left ventriculography. Select a device 0-2mm larger than the VSD size as assessed. The waist diameter Φ D1 of the four small models should be 0-1mm larger than the VSD size, while the waist diameter Φ D1 of the four large models should be 1-2mm larger than the VSD size.
- ◇ Access the VSD with a guidewire using either venous or arterial approach. The venous approach may be from either the inferior vena cava or the superior vena cava depending on the location of the VSD. Form an arteriovenous loop with the wire.
- ◇ Advance the sheath and dilator over the guidewire until the 0.038-inch "J" tip of the dilator crosses the VSD. Then remove the dilator.
- ◇ Advance the guidewire into the apex of the ventricle and position the tip of the sheath in the body of the ventricle.
- ◇ Pass the delivery cable through the loader and attach the occluder device to the tip of the delivery cable by rotating the device clockwise until it is firmly attached. To ensure proper occluder release, rotate the device counter-clockwise 1/8 of a turn.
- ◇ Immerse the occluder and loader in saline solution and retract the device into the loader. Slowly remove the guidewire and allow for back-bleeding to purge air from the system.
- ◇ Connect the loader to the sheath. Transfer the device from the loader into the sheath, and without rotation, advance the occluder to the tip of the sheath.
- ◇ Retract the sheath slowly until the distal disc is deployed. Pull the entire assembly (delivery cable and sheath) into the VSD.
- ◇ Once position is confirmed, retract the sheath to deploy the proximal disc.
- ◇ If the position of the device is unsatisfactory, stabilize the delivery cable and re-advance the sheath until the device is completely within the sheath, and reposition the device and deploy it again, or remove the device from the patient.
- ◇ If the device position is satisfactory, attach the plastic vise to the delivery cable, then release the device by rotating the delivery cable counter-clockwise until it separates from the occluder, finally remove the delivery cable and sheath from the patient.

- ◇ Use TEE to confirm device placement, assess for residual shunts, obstruction, or regurgitation induced by the device. Perform an angiogram to assess for residual flow through the device.

9. Post-procedure

- All patients should be kept overnight for observation.
- Anticoagulation treatment with heparin for 24 hours and antibiotic for 3 days.
- Oral aspirin 3-5 mg/(kg.d) for children, and 3 mg/(kg.d) for adults for 6 months; conditionally, plus plavix 75 mg/d for some adult patients.
- Re-examination by echocardiography, ECG and chest X-ray at 24 hours, 1, 3, 6 and 12 months respectively, after the procedure.

10. Expiry Date

- The KONAR-MF™ VSD Occluder is sterilized by ethylene oxide. The shelf life is indicated on the labels, and the manufacturing date and use by date are marked on the label. DO NOT use an expired product.

11. Package and Label

- The KONAR-MF™ VSD Occluder is supplied sterile.
- The KONAR-MF™ VSD Occluder is fixed and protected on a PETG tray sealed with Tyvek paper, and then the product is sealed with another Tyvek pouch, on which a label and a sterilization indicator are attached. The product is sterilized and put in a box with IFU, implant card, customer feedback form, and certificate of compliance. An outer label is attached on the box.

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Lifetech Scientific (Shenzhen) Co., Ltd
LifeTech Scientific (Europe) B.V.