MagicTouch AVF
SIROLIMUS COATED PTA BALLOON CATHETER

EN	Sirolimus Coated PTA Balloon Catheter (0,035" OTW) Instructions for Use:	Page 04
FR	Sirolimus enduit cathéter à ballon PTA (0,035") Instructions d'utilisation	Page 08
DE	Sirolimus beschichteter PTA Ballonkatheter (0,035") Gebrauchsanweisung	Page 13
IT	Sirolimus rivestito PTA palloncino Catetere (0,035") Istruzioni per l'uso	Page 18
ES	Catéter con globo de PTA recubierto de sirolimus (0,035") Instrucciones de Uso	Page 23
PT	Cateter com balão PTA revestido com Sirolimus (0,035") Instruções para uso	Page 28
RU	БАЛЛОННЫЙ КАТЕТЕР ДЛЯ ЧТА ПОКРЫТЫЙ СИРОЛИМУСОМ (0,035" OTW) Инструкция по использованию	Page 33
TR	Sirolimus Kaplı PTA Balon Kateter (0,035") Kullanım Kılavuzu	Page 38

Symbol used in labeling

Length
Longueur
Länge
Lunghezza
Longitud
Comprimento
длина
Uzunluk

Do not re-use
Ne pas réutiliser
Nicht wiederverwenden
Non riutilizzare
No reutilizar
Não reutilizar
Не использовать повторно
Tekrar kullanmayın

Diameter
Diamètre
Durchmesser
Diametro
Diámetro
Diámetro
диаметр
Çap

Use by Date
Date limite d'utilisation
Herstellungsland
Data di scadenza
fecha de caducidad
Data de validade
срок годности
Son Kullanma Tarihi

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MagicTouch PTA

SIROLIMUS COATED PTA BALLOON CATHETER

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1.0 DEVICE DESCRIPTION

Magic Touch PTA Sirolimus Coated PTA Balloon Catheter (Magic Touch PTA) consists of an over-the-wire (OTW) balloon catheter coated with a polymer free formulation containing Sirolimus drug as an active ingredient in a phospholipid excipient. It uses proprietary Nanolute Coating technology. Generic name of Magic Touch PTA is Sirolimus Coated PTA Balloon Catheter with GMDN Code 62551. Sirolimus Coated PTA Balloon Catheter is single use, minimally invasive, non-implantable drug delivering device intended to use for Percutaneous Transluminal Angioplasty of de-novo, stenotic and in-stent restenotic lesions in superficial femoral, popliteal arteries, below the knee arteries and iliac arteries with reference vessel diameters of 2.00 mm to 12.00 mm.

1.1 PTA CATHETER DESCRIPTION

Sirolimus Coated PTA Balloon Catheter (0.014", 0.018", 0.035" OTW)

Sirolimus Coated PTA Balloon Catheter (0.014", 0.018", 0.035" OTW) is considered a drug/device combination product intended to use for percutaneous transluminal angioplasty of de-novo, stenotic and in-stent restenotic lesions in superficial femoral, popliteal arteries, below the knee arteries and iliac arteries. Sirolimus Coated Balloon Catheter consists of PTA balloon catheter coated with an API (sirolimus) and a phospholipid excipient in a 1:1 ratio.

Sirolimus Coated PTA Balloon Catheter (0.014", 0.018", 0.035" OTW) is Over-the-Wire catheter (OTW) dual lumen catheter with a distally mounted semi-compliant inflatable balloon and an atraumatic tapered tip. The catheter has an OTW construction with a dual lumen proximal shaft. One lumen is used for inflating the balloon and the second lumen allows access to the distal tip of the guide wire insertion. Two marker bands located on the inner member are used to indicate by fluoroscopy balloon's positioning. The tip is soft and tapered to facilitate advancement of the catheter to and through the stenosis. The hub is a Y-arm luer adhesively bonded to the proximal end of the shaft to allow for entry to the guide wire lumen and to allow for connection of the inflation device to the inflation/deflation lumen. The hub is also made with a transparent material to facilitate purge. A strain relief, which is clipped on the luer, is used for device holding and identification. The balloon catheter is covered with a hydrophilic coating except atraumatic tapered distal tip and balloon portion to reduce friction.

The catheter working length are 45, 80, 90, 120, 130, 150 cm and is available in diameters between 1.50 and 12.00mm, and lengths of 20, 40, 60, 80, 100, 120, 150 and 200 mm. The device is provided sterile by ethylene oxide sterilization.


Figure 1: Schematic of Sirolimus Coated PTA Balloon Catheter (0.014"OTW, 0.018"OTW, 0.035"OTW)

(1) PTA Balloon with drug coating; (2) Ring Marker; (3) Distal Tubing; (4) Soft Tip Tubing; (5) Proximal Shaft Tubing; (6) Strain Relief Tubing; (7) Connector; (8) UV Glue; (9) Hydrophilic Coating

Table 4: Specification- Sirolimus Coated PTA Balloon Catheter


Attribute	Sirolimus Coated PTA Balloon Catheter (0.014", 0.018", 0.035" OTW)
Balloon Diameter (mm)	1.50, 2.00, 2.50, 3.00, 3.50, 4.00, 5.00, 5.50, 6.00, 7.00, 8.00, 9.00, 10.00, 12.00
Balloon length (mm)	20, 40, 60, 80, 100, 120, 150, 200
Shaft length (Usage catheter length)	45, 80, 90, 120, 130, 150 cm
Guide wire compatibility	0.014", 0.018", 0.035"
Guiding catheter Compatibility	Minimum of 4F to 5F (depending on balloon size)
Catheter Configuration	Over-the-Wire (OTW)
Tip shape	Tapered
Balloon material	Polyether block amide / PA12
Ring Marker	90% Pt / 10% Ir
Coating on catheter shaft only	Hydrophilic PVP
Protection Wire	Stainless Steel
Balloon Coating	Encapsulated Sirolimus in phospholipid
Drug	Sirolimus
Drug dose	1.27 µg/mm ²
Primary Packaging	Tyvek Pouch
Sterilization	EO
Secondary Packaging	Aluminium Pouch
Box	Cardboard paper

1.2 DRUG COMPONENT DESCRIPTION

The active pharmaceutical ingredient on the Magic Touch PTA is Sirolimus. Sirolimus is a white powder, produced by a Streptomyces Hygroscopius, with the empirical formula C₅₁H₇₉NO₁₃ and a molecular

weight of 914.2g/mole. It is insoluble in water, and melts at approximately 183°C -185°C. The chemical name for Sirolimus is: (3S, 6R, 7E, 9R, 10R, 12R, 14S, 15E, 17E, 19E, 21S, 23S, 26R, 27R, 34aS) -9, 10, 12, 13, 14, 21, 22, 23, 24, 25, 26, 27, 32, 33, 34, 34a -Hexadecahydro -9, 27 - dihydroxy -3 - [(1R) - 2 - [(1S, 3R, 4R) - 4 - hydroxy - 3 - methoxycyclohexyl] - 1 - methylethyl] - 10, 21 - dimethoxy - 6, 8, 12, 14, 20, 26 - hexamethyl - 23, 27 - epoxy - 3 Hpyrido [2, 1 - c] [1, 4] oxazacyclopentriacontine - 1, 5, 11, 28, 29 (4H, 6H, 31H) - pentone; Sirolimus CAS Registry number is 53123-88-9.

Sirolimus has the following chemical structure:

The drug coating of Magic Touch PTA is a non-polymer based formulation, consisting of Sirolimus as the active pharmaceutical ingredient and phospholipid excipients which act as the drug carrier; the coating covers the balloon portion of the Magic Touch PTA. Sirolimus (API) is a macro cyclic lactone produced by Streptomyces Hygroscopius. Sirolimus has potent immunosuppressive and anti-proliferative properties.

The Sirolimus coating is evenly distributed across the working length of the balloon at a surface concentration of 1.27 µg /mm². The key functional characteristic of the formulation is to allow for release of Sirolimus to the tissue of the vascular wall during inflation.

Table: 1 to 3 represents the balloon sizes and the nominal total quantity of Sirolimus on each balloon based on the surface concentration of 1.27 µg/mm².

2.0 INDICATIONS FOR USE

Sirolimus Coated PTA Balloon Catheter is single use, minimally invasive, non-implantable drug delivering device intended to use for Percutaneous Transluminal Angioplasty of de-novo, stenotic and in-stent restenotic lesions in superficial femoral, popliteal arteries, below the knee arteries and iliac arteries with reference vessel diameters of 2.00 mm to 12.00 mm.

3.0 CONTRAINDICATIONS

Use of Magic Touch PTA is contraindicated in the following conditions:

- Patients with a hypersensitivity to Sirolimus drug or its structurally related compounds and/or phospholipid.
- Patients who are not fit for recommended anti-platelet and/or anticoagulant therapy.
- Pregnancy and lactating (breastfeeding) women
- Lesion not suitable for interventional treatment.

4.0 WARNINGS

- Contents supplied are sterilized using Ethylene Oxide (EO) process. Do not use if sterile barrier is damaged or opened prior to intended use. Do not use if packaging or product damage is evident.
- Do not use after the "Use by" date.
- Magic Touch PTA is designed and intended for single use only. DO NOT RESTERILIZE AND/OR REUSE. Reuse or re-sterilization may create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Reuse or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness and death. Concept will not be responsible for any direct, incidental or consequential damages resulting from re-sterilization or reuse.
- Inspect the device prior to procedure, to verify functionality and any sign of damaged parts. Do not use the device if the outer or the inner package is damaged or opened.
- When the Magic Touch PTA is in the body, it should be manipulated only under sufficient and/or high-quality fluoroscopy. Prior to withdrawing the catheter from the lesion, the balloon must be completely deflated under vacuum. If resistance is met during manipulation, determine the cause of the resistance before proceeding.
- Do not use air or any gaseous medium to inflate the balloon. Use only the recommended inflation medium.
- Do not expose the device to organic solvents, e.g. alcohol.
- Do not exceed the Rated Burst Pressure (RBP). Use of a pressure monitoring device is recommended to prevent over pressurization. Inflation in excess of the rated burst pressure may cause the balloon to rupture.
- Magic Touch PTA should only be used by physicians trained in

peripheral vascular percutaneous interventional procedures. Physicians should keep themselves informed and updated on recent publications about PTA techniques.

- To reduce the risk of vessel damage, the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis.
- The safety and effectiveness of Magic Touch PTA has not been established in cerebral, carotid, or renal vasculature.
- Assess risks and benefits before treating patients with a history of uncontrollable allergies to contrast agents.
- Use of more than two Magic Touch PTA deployed in a single target lesion during a single procedure has not been assessed.
- Handle the product with caution to avoid any damage to the balloon coating or folded balloon.
- Take appropriate precautions to prevent or reduce clotting when catheter is used. Flush or rinse all products (except Magic Touch PTA Balloon) entering the vascular system with sterile isotonic saline or a similar solution via the guide wire access port prior to use. Consider the use of systemic heparinization.
- Prior to the procedure, inspect the product to verify that the product is intact. DO NOT use a catheter that has been damaged.
- To minimize the possible introduction of air into the system, it's imperative that, before proceeding, careful attention is paid to the maintenance of tight catheter connections and through aspiration and flushing of the system.
- Magic Touch PTA should be used with caution for procedures involving calcified lesions due to the abrasive nature of these lesions.
- Never advance the angioplasty catheter without the guide wire extending from the tip.
- Never attempt to move the guide wire when the balloon is inflated. Do not advance the catheter against significant resistance. The cause of resistance should be determined via fluoroscopy and appropriate remedial action taken.
- This product is not intended for the expansion or delivery of a stent.
- The use of this product carries the risks associated with percutaneous transluminal angioplasty, including thrombosis, vascular complications, and/or bleeding events.
- As on date there are neither long term clinical studies performed, nor well controlled studies been published to evaluate the carcinogenic or genotoxic potential of the drug sirolimus coated balloons. Therefore, it is uncertain at present whether API Sirolimus (Rapamune, Rapamycin) is safe to be used in any kind of treatment during pregnancy or in patients intending to have children.

5.0 PRECAUTIONS

5.1 Uses in Conjunction with Other Procedures

- It is not recommended that the Magic Touch PTA be used in conjunction with other drug coated balloons or drug eluting stents to treat the same lesion in the same procedure. The safety of Magic Touch PTA in conjunction with different drug-device products has not been assessed.
- Use of the Magic Touch PTA has not been studied in conjunction with other interventional techniques such as atherectomy, laser, or cutting/scoring balloons.
- Potential interactions of Magic Touch PTA with alternative therapies such as drug-eluting stents, lasers, atherectomy, cryoplasty, cutting/scoring balloons, and brachytherapy has not been evaluated and should be avoided whenever possible.

5.2 Uses in Special Population

Paediatric Use: The safety and effectiveness of Magic Touch PTA in paediatric patients (under 18 years of age) has not been established.

5.3 Device Handling Precautions

- Do not immerse Magic Touch PTA in a saline bath. Avoid any pre-insertion contact with fluid on the balloon as it may affect the integrity of the drug coating. Replace any Device where the balloon has come into contact with fluids prior to use.
- The coated balloon portion of Magic Touch PTA must be handled with dry sterile gloves. Care should be taken to minimize unnecessary contact with the coated balloon portion of Magic Touch PTA during preparation and insertion.
- The balloon protector and wire lumen stylet should stay in place during preparation of the Magic Touch PTA and not be removed until just prior to placing over guide wire.
- If difficulty is encountered while removing the balloon protector after flexing, a new Magic Touch PTA should be utilized.

5.4 Device Use/Procedure Precautions

- To ensure therapeutic drug delivery:
 - Never inflate Magic Touch PTA prior to reaching the target lesion.
 - Magic Touch PTA should be advanced to the target site in an efficient manner (≤ 3 minutes) and immediately inflated.
 - Maintain balloon inflation for a minimum of 60 seconds however, inflation time up to 180 seconds is highly recommended to get optimum drug transfer and results. Use the maximum balloon inflation time as per your institution's standard of care.
- After insertion, do not over-tighten the haemostatic adaptor (if used) around Magic Touch PTA shaft as lumen constriction may occur, affecting inflation/deflation of the balloon.
- Always advance and retrieve Magic Touch PTA under negative pressure.
- The Magic Touch PTA should always be manipulated under fluoroscopic observation when in the body.

- Do not continue to use Magic Touch PTA if the shaft is bent or kinked.
- Whenever possible, the Magic Touch PTA should be the final treatment of the vessel; however, post-dilatation is allowed with other PTA catheter or used Magic Touch PTA.

5.5 Pre- and Post-Procedure Antiplatelet Regimen

Dual antiplatelet therapy should be administered according to current standards pre-procedure and for a minimum of 4 weeks after the intervention. Prolonged antiplatelet therapy can be given at the discretion of the physician.

6.0 DRUG INFORMATION

6.1 Mechanism of Action

Magic Touch PTA coating contains Sirolimus, an anti-proliferative pharmaceutical agent that specifically binds to and stabilizes microtubules. Sirolimus is a powerful inhibitor of T-Lymphocyte activation, smooth muscle cell and endothelial cell proliferation in response to stimulation by cytokines, growth factors and antigens. In cells, sirolimus binds to the immunophilin, FK Binding protein-12 (FKBP-12). The sirolimus- FKBP- 12 complex binds to and inhibits the activation of the mammalian Target of Rapamycin (mTOR), a key regulatory kinase. This inhibition suppresses cytokine-driven cell proliferation, inhibiting the progress of the cell cycle from the G₁ to the S phase.

6.2 Drug Interaction

- Drug interaction studies have not been performed with Magic Touch PTA because of limited systemic exposure to Sirolimus eluted from Magic Touch PTA. However, consideration should be given to the potential for both local and systemic drug interactions in the vessel wall when deciding to implant Magic Touch PTA in a subject taking a drug with the known interaction with Sirolimus or when deciding to initiate therapy with such a drug in a patient who has recently undergone a procedure with Magic Touch PTA.
- Sirolimus is metabolized by cytochrome P450 3A4 (CYP3A4) in the gut wall and liver and undergoes efflux from enterocytes of the small intestine by p-glycoprotein (P-gp). Absorption and subsequent elimination of systemically absorbed sirolimus may be influenced by drugs that affect CYP3A4 and P-gp.
- Medications that are strong inhibitors of CYP3A4 and P-gp may increase Sirolimus levels while inducers of CYP3A4 and P-gp may reduce Sirolimus metabolism in vivo.
- Sirolimus when prescribed as oral medication may interact with the drugs/foods listed below:
 - P-gp Inhibitors such as but not limited to Cyclosporine, Digoxin
 - Antibiotics such as but not limited to Ciprofloxacin, Ofloxacin, Rifabutin
 - Glucocorticoids
 - CYP3A4 Inhibitors such as but not limited to Ketoconazole, Fluconazole
 - CYP3A4 Inducers such as but not limited to Rifampicin, Phenytoin, Dexamethasone
 - Herbal preparations (St. John's wort- Hypericum Perforatum)
 - Amphotericin B such as but not limited to Abcetyl, Amphotin
 - Cimetidine such as but not limited to Tagamet
 - Vitamins
 - Vaccination- Sirolimus is potent immunosuppressant hence it may affect response to vaccination making the vaccination less effective. For some period after receiving Magic Touch PTA, use of live vaccines such as, but not limited to measles, mumps, BCG, yellow fever, varicella, Ty21a typhoid should be avoided.
- Grapefruit/grapefruit juice reduces CYP3A4 mediated metabolism of Sirolimus.
- Medicines to lower cholesterol or triglycerides
- Medicines for high blood pressure or heart problems
- Anti-seizure medicines
- Medicines used to treat stomach acid, ulcers, or other gastrointestinal problems
- Drugs that may increase Sirolimus blood concentration:
- Antifungal agents such as but not limited to Itraconazole, Cotrimoxazole
- Gastrointestinal prokinetic agents such as but not limited to metoclopramide
- Macrolide antibiotics such as but not limited to Troleandomycin
- Calcium channel blockers such as but not limited to Diltiazem, verapamil
- Other drugs like Danazole, HIV Protease Inhibitor (Indinavir)
- Drugs that may decrease sirolimus level
- Anticonvulsants such as but not limited to Phenobarbital

7.0 POTENTIAL ADVERSE EVENTS

Potential adverse events which may be associated with a peripheral balloon dilatation procedure include additional intervention, Allergic reaction to drugs or contrast-medium, Amputation/loss of limb, Aneurysm or pseudo aneurysm, Arrhythmias, Arterial spasm, Arteriovenous fistula, Vessel dissection, perforation, rupture, or spasm, Embolization, Haematoma, Haemorrhage, including bleeding at the puncture site, Histological changes in vessel wall, including inflammation, cellular damage or necrosis, Hypotension/hypertension, Nausea & vomiting , Occlusion, Pain or tenderness, Pneumothorax or hemothorax, Sepsis/infection, Shock, Stroke, Thrombosis, Vascular complications which may require vessel repair. Although systemic effects are not anticipated, refer to the Physicians' Desk Reference for more information on the potential

adverse events observed with Sirolimus. There may be other potential adverse events that are unforeseen at this time.

8.0 PATIENT COUNSELING INFORMATION

Physicians should consider the following in counselling patients about this product:

- Discuss the risks associated with peripheral vascular percutaneous interventional procedures.
 - Discuss the risks associated with a Magic Touch PTA
 - Discuss the risks/benefits issues for the particular patient.
 - Post-procedure antithrombotic regimen.
 - Discuss alteration to current lifestyle immediately following the procedure and over the long term.
- #### 9.0 HOW SUPPLIED
- **STERILE:** This device is Non-pyrogenic and sterilized with Ethylene Oxide gas.
 - **WARNING: Do not re-sterilize and/or reuse or reprocess the Device. Do not use if sterile barrier is damaged or opened prior to intended use. This Device is for single patient use only.**
 - **CONTENTS:** One (1) Magic Touch PTA [Over-The-Wire (OTW)].
 - **STORAGE:** Store in a dry, dark place. Store between 8°C to 25°C and use before the "Expiry date". Do not store near radiation or ultra-violet light sources.
 - **DISPOSAL INSTRUCTIONS:** After use, this product may be a potential biohazard. Handle and dispose-off in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.
- #### 10.0 DIRECTIONS FOR USE
- #### 10.1 Materials Required
- In addition to Magic Touch PTA, the following standard materials may also be required:
- Introducer sheath
 - Guide catheter
 - Contrast-medium
 - Sterile saline
 - Inflation device with manometer
 - Leuer-lock syringe for purging
- #### 10.2 Inspection Prior to Use
- Prior to angioplasty, carefully examine all equipment to be used during the procedure, including sterile package and dilatation catheter to verify proper function. Do not use after the "Use By" date. Verify that Magic Touch PTA and sterile packaging have not been damaged in shipment transit.
- Warning: Contents supplied are sterilized using Ethylene Oxide (EO) process. Do not use if sterile barrier is damaged or opened prior to intended use.**

- #### 10.3 Use of Multiple Magic Touch PTA Sirolimus Coated PTA Balloon Catheters
- The extent of the patient's exposure to the drug coating is directly related to the number of balloons used. Refer Size Matrix and Sirolimus Content (Table 1 to 3) for Magic Touch PTA containing the nominal Sirolimus content for each device size.
- If multiple Magic Touch PTA are required due to a lesion length greater than the longest available Magic Touch PTA, the balloons must overlap by at least 2.5 mm. The size of additional Magic Touch PTA should not be longer than required to allow for this overlap and complete the lesion coverage with about 2.5 mm extended beyond the lesion both proximally and distally. Proper size selection is important to avoid excessive overlap.
- Warning:** The safety and effectiveness of implanting multiple Magic Touch PTA in a patient has not been clinically evaluated.
- #### 10.4 Magic Touch PTA Sirolimus Coated PTA Balloon Catheter Preparation
- Important note: Pre-dilatation is necessary before use of Magic Touch PTA.

1. Remove the device from the packaging.
2. Verify the balloon size is suitable for the procedure and the selected accessories are compatible with the catheter as labelled.
3. Guide Wire Lumen Flush:
Flush the guide wire lumen with HepNS until fluid exist the distal end of the delivery system.
4. Prepare the inflation device/syringe with diluted contrast medium.
Warning: Use the recommended balloon inflation medium of contrast and sterile saline (≤ 50% contrast). Never use air or any gaseous medium to inflate the balloon.
5. Prior to use, the air in the balloon catheter should be removed. To facilitate purging, select a syringe or inflation device with a 10 ml or larger capacity and fill approximately half of it with the recommended contrast medium.
6. Connect a stopcock to the balloon inflation Leuer hub on the dilatation catheter.
7. Connect the syringe to the stopcock.
8. Hold the syringe with the nozzle pointing downward, open the stopcock and aspirate for approximately 15 seconds. Release the plunger.
9. Repeat from step 4 above as needed until bubbles no longer appear during aspiration (negative pressure). Once completed, evacuate all air from the barrel of the syringe/inflation device.

10.5 Use of Magic Touch PTA Sirolimus Coated PTA Balloon Catheter

1. While under negative pressure and before removing the balloon protecting cap, perform the following steps to reduce friction between the balloon protecting cap and the balloon and remove the balloon protector:
 - Step 1- Grasp the proximal end of the balloon protecting cap with one hand.
 - Step 2- Using the opposite hand, gently slide the thumb and forefinger from the proximal end of the balloon protecting cap out toward the distal end of the balloon protecting cap while flexing the balloon slightly downward approximately 15 degrees.
 - Step 3- Repeat Step 2, but flex the balloon protecting cap upward approximately 15 degrees.
 - Step 4- Grasp the balloon protecting cap at roughly the midpoint of the balloon protecting cap and pull it away from the balloon catheter.
2. With the catheter tip oriented down/vertically, flush the wire lumen.
3. Backload the distal tip of the dilatation catheter onto the guidewire.
4. While the balloon is still fully deflated and under negative pressure, slowly advance the Magic Touch PTA through the introducer sheath and over-the-wire to the site of inflation. During catheter advancement, inspect the catheter shaft for damage.
5. To ensure therapeutic drug delivery, the Magic Touch PTA should be advanced to the target site in an efficient manner (≤ 3 minutes) and immediately inflated.
Warning: Do not exceed the Rated Burst Pressure (RBP) recommended for this device. Balloon rupture may occur if the RBP rating is exceeded. To prevent over-pressurization, use of a pressure monitoring device is recommended.
6. Position the balloon relative to the lesion, ensuring coverage of at least 5mm proximally and distally beyond the margins of the lesion segment, and inflate the balloon to the appropriate pressure (reference Compliance Chart). The use of a radiopaque ruler is recommended to ensure appropriate placement of the Magic Touch PTA.
7. Apply negative pressure to fully deflate Magic Touch PTA. Prior to removal, confirm that the balloon is fully deflated under fluoroscopy.
8. Perform angiography to confirm dilatation of the lesion.
9. Withdraw Magic Touch PTA from the body under negative pressure. Maintain the guidewire across the stenosis.
10. Whenever possible, Magic Touch PTA should be the final treatment of the vessel; however, postdilatation is allowed with another PTA catheter or used Magic Touch PTA.
11. After confirming that a satisfactory dilatation was achieved, remove all equipment from the body and close access site per standard clinical practice.
12. Refer Pre- and Post-Procedure Antiplatelet Regimen for the dual antiplatelet pharmacological therapy recommended with use of Magic Touch PTA.
13. After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and local applicable laws and regulations.

11.0 PATIENT INFORMATION

In addition to these Instructions for Use booklet, a Product treatment Card is provided in the package. It is advised that patients keep this card in their possession for procedure / Balloon Catheter identification.

12.0 DISCLAIMER OF WARRANTY

Concept Medical B.V. has endeavoured to ensure that the products comply with all relevant standards and regulations currently in force and to ensure that the quality of the products meets the requirements of the above-mentioned standards and regulations for a period ending upon the indicated expiry date. The above statement does not apply where the products are used for a purpose other than its intended purpose. Where any loss or damage is caused (other than death or personal injury) due to a defective product, Concept Medical B.V. shall not be liable for such loss or damage.

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