0344

**EN****AcuMark™ Sizing Balloon  
Instructions for Use****BG****Балон за оразмеряване AcuMark™  
Инструкции за употреба****DA****AcuMark™ dimensioneringsballon  
Brugsanvisning****EL****Καθετήρας με μπαλόνι προσδιορισμού  
μεγέθους AcuMark™  
Οδηγίες χρήσης****FI****AcuMark™-pallokatetri defektin mittaamiseen  
Käyttöohjeet****HR****Balon za dimenzioniranje AcuMark™  
Upute za uporabu****IT****Palloncino di dimensionamento AcuMark™  
Istruzioni per l'uso****NL****De AcuMark™ maatbepalingsballon  
Gebruiksaanwijzing****PL****Balon do wymiarowania AcuMark™  
Instrukcja obsługi****RO****Balonul de dimensionare AcuMark™  
Instrucțiuni de utilizare****SV****AcuMark™ dimensioneringsballong  
Bruksanvisning****AR****AcuMark™ بالون التحجيم  
تعليمات الاستخدام****CS****Balónkový katétr pro měření AcuMark™  
Návod k použití****DE****The AcuMark™ Größenbestimmungsballon  
Gebrauchsanweisung****ES****Catéter-balón de medición AcuMark™  
Instrucciones de uso****FR****Ballonnet de calibre AcuMark™  
Mode d'emploi****HU****Az AcuMark™ méretezőballon  
Használati utasítás****LT****„AcuMark™“ dydžio nustatymo balionas  
Naudojimo instrukcijos****NO****AcuMark™ dimensjoneringsballong  
Bruksanvisning****PT****Balão de dimensionamento AcuMark™  
Instruções de utilização****RU****Измерительный баллонный катетер  
AcuMark™  
Инструкция по применению****TR****AcuMark™ Ölçme Balonu  
Kullanım Talimatları****Lifetech Scientific (Shenzhen) Co., Ltd.**

# Content

EN: Instructions for Use .....	1
AR: تعليمات الاستخدام .....	6
BG: Инструкции за употреба.....	10
CS: Návod k použití.....	16
DA: Brugsanvisning.....	21
DE: Gebrauchsanweisung.....	26
EL: Οδηγίες χρήσης .....	32
ES: Instrucciones de uso.....	38
FI: Käyttöohjeet .....	43
FR : Mode d'emploi .....	48
HR: Upute za uporabu.....	54
HU: Használati utasítás.....	59
IT: Istruzioni per l'uso.....	64
LT: Naudojimo instrukcijos.....	70
NL: Gebruiksaanwijzing.....	75
NO: Bruksanvisning.....	80
PL: Instrukcja obsługi .....	85
PT: Instruções de utilização.....	90
RO: Instrucțiuni de utilizare .....	95
RU: Инструкция по применению.....	100
SV: Bruksanvisning .....	106
TR: Kullanım Talimatları .....	111

## 1 Device Description

The AcuMark™ Sizing Balloon is a triple lumen catheter which is designed to measure cardiovascular defects (such as an ASD). The balloon is made of soft compliant TPU (Thermoplastic Polyurethane) plastic membrane. Three radiopaque marker bands located inside the balloon at 5 mm (as measured from inner edge to inner edge between adjacent marker bands) intervals are used as a distance reference for magnification correction (Figure 1).


Figure 1 Illustration scheme of AcuMark™ Sizing Balloon

## 2 Device Specification and Size Parameters

Table 1 The specification of the AcuMark™ Sizing Balloon

Catalogue NO.	Balloon Length(m m)	Maximum Defect Size(mm)	Shaft Size	Shaft Usable Length (mm)	Guidewire(in)	Maximum Inflation Volume(cc)
LT-SZB-24	40	27	7F	750	0.035	30
LT-SZB-34	50	42	8F	750	0.035	90

Table 2 The compliance of the AcuMark™ Sizing Balloon

Catalogue NO.	Inflation volume(cc)		Balloon Diameter (mm)
	Nominal	Maximum	
LT-SZB-24	Nominal	19	24
	Maximum	30	29
LT-SZB-34	Nominal	47	34
	Maximum	90	44

## 3 Indications for Use

The AcuMark™ Sizing Balloon is intended for use in those patients with cardiovascular defects wherein accurate measurements of the defect is important to select an appropriately sized occlusion device. Measuring can be accomplished in two ways: radiographic or echocardiography.

## 4 Contraindications

- Patients with hemorrhagic physical disease
- Patients unable to receive or limited to antiplatelet and/or anticoagulation therapy

- Pregnant or breastfeeding women, and women planning to become pregnant
- Known hypercoagulable state or any coagulation disorder
- Allergic or intolerant to contrast media or drugs used for intravascular surgery

## 5 Precautions and Warnings

Read the instructions for use carefully prior to using the device.

- The device should only be used by those physicians trained and experienced in transcatheter defect sizing techniques.
- The AcuMark™ sizing balloon is intended to be used in an operating room or a hospital catheterization lab or clinics, by professionally trained clinicians, such as cardiologists or cardiovascular surgeons. There are unpredictable risks associated with its use in environments other than those described above.
- Do not insert the sizing balloon through an introducer sheath.
- This device is for single use only. Do not reuse. Reuse may cause patient infection or cross-infection.
- Do not use this device if the packaging has been opened or damaged or the sterile barrier is compromised in any way.
- Do not exceed the recommended maximum volume of inflation for each size (see Table 1).
- During the procedure, patient should be heparinized, the device surface shall be flushed with sterile heparinized saline.
- Any patient known to have a bleeding disorder, untreated ulcer or any other contraindications to aspirin therapy, unless another anti-platelet agent can be administered for a minimum of 6 months.
- This device shall be used before the expiry date.
- Thrombus caused by air should be avoided when delivering the sizing balloon to the target position.
- Immediately replace new sizing balloon when the shaft is bent or damaged during the procedure.
- This balloon catheter is not recommended for pressure measurement or fluid injection.

## 6 Potential Complications

The complications that may result from a balloon procedure include:

- Vessel perforation/dissection
- Thromboembolism
- Infection
- Vessel spasm
- Hematoma
- Embolization or Vascular thrombosis
- Pain and tenderness
- Endocarditis
- Arrhythmias
- Short-term hemodynamic compromise
- Valvular interference
- Expansion of the defect

- Cardiac obstruction

Potential complication specific to balloon sizing include:

- Movement of the balloon toward the mitral valve or right atrium
- Enlargement of the ASD
- Obstruction of venous return from the Inferior Vena Cava

**Note:**

- For any serious incident, the user or patient should report it to the manufacturer and the local competent authority (CA).
- The patient should contact his/her doctor whenever they have any questions, concerns or feel any discomfort.

## **7 Directions for Use**


### 7.1. Pre-procedure

- ◆ Patients should be fully heparinized throughout the procedure with a minimum Active Clotting Time (ACT) of 200 seconds prior to device insertion.
- ◆ Following percutaneous puncture of the femoral vein, perform a standard right heart catheterization, as per hospital protocol.
- ◆ Perform an angiogram in order to demonstrate the atrial communication. Catheterize the left atrium through a 45°LAO position and a 35-45°cranial angulation, inject contrast medium into the left atrium or into the right upper lobe pulmonary vein.

### 7.2 Operation


- ◆ Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken. Remove the device from the packaging and remove the balloon cover carefully to avoid any sharp bends or kinks to the balloon.
- ◆ Inflate the balloon with normal saline and exhaust air as shown in Figure 2, or inflate and deflate with carbon dioxide. Flush the other cavity with heparinized saline. Check inflation /deflation time, repeat the procedure several times to verify the inflation/deflation time.

**Note: Refer to table 2 for the range of inflation volume of the balloon.**

**Figure 2.**

- ◆ Procedures will be carried out according to the type of defects or measurements. As the balloon is a soft, compliant TPU (Thermoplastic Polyurethane) plastic membrane, there is no need for expansion of the puncture site. To facilitate this percutaneous entry, forceful negative pressure should be applied to the balloon catheter with an attached syringe.
- ◆ Enter the vessel percutaneously using the standard Seldinger technique over an 0.035" exchange "J" tip stiff guidewire (can be assisted by the use of a multipurpose catheter) positioned in the left upper pulmonary vein or left atrium.
- ◆ Under fluoroscopic and echocardiographic guidance, the balloon catheter will be placed across the defect until the middle radiopaque marker band of the three radiopaque marker bands passes through the defect, as shown in Figure 3. The balloon will then be filled with diluted contrast agent (3:1 ratio) using indeflator until color flow shunt ceases as observed by echocardiography and a waist appears in the balloon.

**Figure 3**

- ◆ Using echocardiography, adjust the image for best visualization, and measure the minimum diameter between the septal wall contact points with the balloon. If a balloon waist is detected, radiographic measurement may also be used to confirm defect diameter.
- ◆ The diameter of this waist is measured by TEE or angiogram. This size corresponds to the stretched size of the defect.

**Note:** Angiogram measurement needs to be corrected for the magnification. The magnification correction can be checked by using the radiopaque marker bands on the catheter shaft.

- ◆ Do not remove the guidewire from the sizing balloon at any time during the procedure to avoid kinking and/or migration.
- ◆ Using a syringe, deflate the balloon completely for 90s.
- ◆ Gently withdraw the sizing balloon from the patient. As the balloon exits the vessel, use a smooth, gentle and steady motion. If you feel resistance when withdrawing the balloon catheter, especially if/when it is known or suspected that the balloon catheter is ruptured or leaked, the balloon catheter should not be forcibly withdrawn. Negative pressure can be applied again for 90s to ensure that the balloon is completely deflated. Inspect the catheter for kinks, then if no kinks or obstructions are detected, withdraw the balloon catheter along with the guidewire under the guidance of fluoroscopic imaging.
- ◆ Apply pressure to the insertion site according to standard practice or hospital protocol for percutaneous vascular procedures.

## 8 Expiry Date

The AcuMark™ Sizing Balloon is sterilized by ethylene oxide. The production date, lot number and expiry date are marked on the label.

**Warning:** DO NOT use expired product.

## 9 Storage and Disposal

- ◆ Keep away from sunlight, keep dry, store between 10 °C-30 °C.
- ◆ Be careful in transportation. Do not squeeze, expose to liquids, bend or expose to extreme heat during transportation.
- ◆ If there are any questions before using the AcuMark™ Sizing Balloon, contact the agent or the manufacturer directly.
- ◆ After use, dispose of the device and packaging in accordance with hospital, administrative and/or local government policy.

## 10 Package and Labeling

The AcuMark™ Sizing Balloon is supplied sterile. The device is immobilized and protected on a scale board, then sealed in double Tyvek1073B pouches, where a primary label is applied. The product is sterilized and inserted into a box with IFU, customer feedback form and certificate of compliance. The outer label is applied on the box. On the primary and outer label, production model, sterilization date, expiry date, lot number, bar code, and symbol for sterilization are clearly indicated.

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LifeTech Scientific(Europe)B.V.